

Perfect medical workers protection

Auto door, foot switch

The automatic visual glass door is adopted for observing the operation inside the tank. The door can be opened by gently kicking the foot switch.

Totally-enclosed

Rider wash-disinfection tank adopts totally-enclosed structure to avoid sanitizer odor leakage so that protects the medical workers.

Automatic adjuvant injection

Medical workers only need to link the adjuvant container with the tank through a connector to avoid harm that caused from manual adjuvant injection.













QPQ 50 series

Parameters		QPQ 50A	QPQ 50B	
Weight (Kg)		100 kg	110kg	
Water consumption/cycle		25L	20L	
Power supply		220V 50Hz		
Power (Kw)		4	4.5	
Water	Water pressure Mpa	0.35~0.4		
	Water hardness	CaCo3 concentration<120mg/L		
	Water PH value	6.5~9		
Water drainage		Upper drainage		

QPQ 60 series

Parameters		QPQ 60A	QPQ 60B		
Weight (Kg)		200kg	220kg		
Water consumption/cycle		25L	20L		
Power supply		220V 50Hz			
Power (Kw)		8	9		
Water	Water pressure Mpa	0.35~0.4			
	Water hardness	CaCo3 concentration<120mg/L			
	Water PH value	6.5~9			
Water drainage		Upper drainage			

万容机械有限公司 WANROOE MACHINERY CO., LTD.

江苏省张家港市人民东路19号国泰新天地广场1311,1313室 Add:Room1311,1313, Guotai New World Plazza, No19, East Renmin Road, Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663
Fax:+ 86 512 58810663
Email:info@wanrooemed.com
Http://wanrooemed.com



Rider Endoscope Washer-Disinfector

For health Be credible





万容机械有限公司 WANROOE MACHINERY CO., LTD.



Rider endoscope washer-disinfector is professional automatic disinfection equipment that developed according to the standard ISO 15883-4.

It can work with various kinds of sanitizer such as: glutaraldehyde, peracetic acid, phthalic dicarboxaldehyde, chlorine dioxide, acidic electrolyzed oxidizing water.







Endoscope protection design

Leaking test function


The leaking test will finish before the endoscope contacts the liquid and will also keep testing during the operation. The system gives out visual and audible signals and will terminate operation as soon as the leaking value exceeds the setting value.

Blocking test function

The blocking test will start before the endoscope contacts the liquid and after the operation. The system gives out visual and audible signals and will terminate operation as soon as the pipe is blocked.

Self-disinfection function

The self-disinfection program can be started during maintenance or after the termination of operation. This function is design to disinfection for wash-disinfection tank (includes 0.2µm filter) and pipeline.



Most speedy wash and disinfection

An endoscope can be wash-disinfected within 15mins by Rider which dramatically improves the endoscope turnover efficiency.



QPQ60 Series

Process traceability system

Process data printing

For filing purpose, the printer will print all of process data timely during the operation.

Process data management

The system collects information of endoscope and operator and connects the process data with computer managing system to achieve synchronization management of patient and disinfection information.





operator information read

endoscope information read

100% wash-disinfection

Omnibearing wash-disinfection

The tank equipped with rotating nozzle and spray arm to wash the outer surface of the endoscope, in the meanwhile, the circulating water flow keeps flushing the endoscope from inside.

Endoscope pressure booster pump

It is used for washing of tiny pipe, such as air/water injection nozzle, to reduce the formation probability of bacteria

Sterile water rinse

To avoid secondary contamination caused by unsanitary water, the endoscope will be rinsed by sterile water which filter by 0.2µm filter after disinfection.

Drying function

It is used for drying the inner and outer surface of endoscope. Rider equipped with air drying and alcohol drying for user to select.



