

BAE500 Surgical bed



Feature:

- StayPut Retractable bed frame design(Patient will feel very convenient while be on drip or used the heart monitor since it won' t need to change patient position when the bed back section rise and low.)-optional
- Central braking system-6 inch diameter castor/fifth castor for direction all reached by foot pedal located in middle of bed.
- Made of superior cold-rolled steel
- One step forming sleeping bed
- Treated by shot blasting, electrophoresis and powder coating
- Four pieces of PP complete side rails, full protection with zero gap, tuck-away
- Nurse controls on footboard and siderails
- Fixed patient controls on siderails
- Electronic function lockout controls on footboard
- Degree indicator for head elevation and leg section
- Four independent motors, standard: LINAK brand
- Automatic Weighting-optional
- CPR quick release handles on both side of bed
- battery back up
- auto contour feature simultaneously raises both the head and the knees
- Wire hidden



Details about the ICU Surgical Bed



In-bed Scale, Accurate, repeatable weighting with the patient any position



Central brake/fifth castor system all in one



CPR quick release handles



Wire hidden



Angle indications



Metal bracket, long use time



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanroomed.com www.wanroomed.com

BAE500 Surgical bed



General Description:

This type of Surgical bed consists of Cold-rolled steel plate and four pcs of long PP siderails. It is mainly used in ICU room. The most advantage is the embedded controller on the foot board and both two sides of the siderails, which allow caregivers to spend their time helping patients recover quickly in a secure and comforting environment.

Standard Specifications

Model number	BAE500
Measurement	Overall size(L*W*H): 2230*1050*450-790mm(with 6 inch diameter castor)
Patient	Touching surface: 1960*900mm
Weight Load	Fully tested sturdy construction capable of taking a maximum user weight of up to 250kg
Function	Backrest max upward angle: 75°
	Footrest max upward angle: 40°
	Height: 450-790mm (castor diameter:6 inch)
	Trend/Reverse Trend. $\pm 12^\circ$
Power	220VAC or 110VAC Frequency: 50-60HZ Power Plug: EU type standard if other type, confirm each order

Standard Accessories List

Headboard/Footboard	1 set
Siderails	4 pcs
I.V pole	1 piece
Castor	4 pieces/5 inch
Drainage Hook	2 pieces/each side, total of 4pcs of hooks
Motor	4 pcs

Optional:

1. Over Bed table
2. beside locker
3. Mattress(sponge mattress, or Bedsore prevention mattress)



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com