

BIC800

Bariatric Care ICU Bed



Safety and Efficiency, Mobility and Dignity

The BIC800 Electric Hospital Bed provided by WANROOEMED is an American style electric adjustable bed that is with good quality, excellent function and affordable price comparing with American made hospital bed. This new adjustable bed for overweight patients provides progressive mobility with flexible wound and pulmonary therapies.

Optional Air Mattress:

Helps prevent and treat pressure ulcers.

Feature:

1. Length extension.
2. One button chair position.
3. Bed frame and leg are all made of cold-rolled steel plate and tube and coated by electrostatic spray after twice phosphatization.
4. Back section is made of American electro plate, X-ray available.
5. Head and foot board, side Rails are all made of imported PP materials. Side rails with angle indication and operation buttons.
6. Four bumpers on easy corner.
7. It is equipped with LINAK motor imported from Denmark which is quiet and noiseless and able to support over 250kgs.
8. Backrest and footrest can move at the same time only by one button on the panel.
9. Centrally controlled silent castors of high strength and high wear resistance enable the bed to move flexibly, lightly and conveniently.
10. CPR release mechanism.
11. Backup battery.
12. Direct wheels, easy to control direction.(optional)
13. With infusion stand and drainage hook.
14. With night light system.
15. Weighting scale system(optional)



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com

BIC800

Bariatric Care ICU Bed



Button control:
Press the button to get what position you want.
Head Angle indication
Monitor patient head angle in any siderail position at point of care

Length extension:
Foot section retractable length allows the bed system to be customized to the patient's height, minimizing sliding down in bed



Battery backup:
No need to find an outlet when repositioning the bed during patient transfer.

Safer brakes—Central lock castor, step the brakes to red color direction to lock the castors, negative direction for release the castor
Fifth castor: with fifth castor to steer the direction easily(optional)



Head section and back section angle indication:
Monitor patient head angle in any siderail position at point of care

Bumper rollers on 4 corners:



CPR release:
A:CPR buttons on the siderails
B:CPR manual handle under the head section

Head section,X-ray available:
Under head section with cassette track for photograph use



Safer weighing and exit monitoring(optional)
Easily weigh your patient in place with superior accuracy - no need to transport to a scale.
Highly accurate, sensitive system. Based on patient weight for fewer false alarms; deploys if the system is compromised.

Good export wooden package, no need fumigation. Film bag pack—carton pack—wooden case. Whole set package, no need complicated installation.



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com