

WANROOEMED

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

AOT100 Electro-hydraulic Operating Table

Developments made in the conjunction with doctors: this is one of the basic principles which has made WANROOEMED the world market leader for operating tables. Close cooperation and years of experience have resulted in tried and tested products which facilitate the work in the operating room and help to increase the safety and comfort of the patient.



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

WANROOEMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

AOT100

Electro-hydraulic Operating Table



江苏如翼医疗科技有限公司
JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.
Add: Room1311,1313, Guotai New World Plaza, No19,
East Renmin Road,Zhangjiagang, Jiangsu, China.
Tel:+86 512 58810663
Whatsapp/Mobile:+ 86-18962287398;15962399164
Email:info@wanrooemed.com
Http://wanrooemed.com

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

WANROOEMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

AOT100

Electro-hydraulic Operating Table



Product Description

AOT100 as an advanced automatic operating table, it can meet various requirements of operations. Tabletop is X-ray available. Extremely low tabletop, especially suitable for cerebral surgery. Engineered for increasingly complex surgical scenes, AOT100 comes with intelligent automatic control system satisfying various demands of general surgeries. Detachable head plate ensure the table can be use with the head stock.

- 520mm minimum and 820mm maximum heights for comfortable working, both seating and standing.
- Memory pad with antistatic, waterproof design.
- Optional battery power for 50-80 operations.



Details

- 1, Leg plate: Imported gas spring for easy controls.
- 2, Two joint head plate.
- 3, Built-in kidney bridge.
- 4, Extreme minimum height 520mm.



TABLE HIGHEST HEIGHT



TABLE LOWEST HEIGHT



WANROOEMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

Model	Site	Length	Width	Tabletop lowest height	Elevation range
AOT100		≥2085mm	≥550mm	≤520mm	300mm

PARAMETER AOT100 electro-hydraulic operating table

Movements		Model: AOT100
Tabletop trendelenburg	Reverse Trendelenburg	≥25°
	Trendelenburg	≥25°
Lateral Tilt (Left And Right)		≥15°
Head plate	Up	≥40°
	Down	≥90°
Back plate	Up	≥80°
	Down	≥20°
Leg plate	Up	≥15°
	Down	≥90°
	Outward	≥90°
Kidney Bridge Elevation		110 mm

STANDARD ACCESSORIES

Anaesthesia screen	1 piece	Body support	1 pair	Leg support	1 pair	Hand remote	1 piece
Kidney bridge handle	1 piece	Arm support	1 pair	Foot support	1 pair	Mattress	1 set
Shoulder support	1 pair	Long fixing clamp	1 pair	Fixing clamp	8 pieces	Power Line	1 piece

