

First Aid Kits



WANROOEMED

FIRST AID KITS

WANROOEMED | YOUR RELIABLE PARTNER

EFK-07C General Aluminum First-aid Kit

Size: 55.7*39*20cm

The EFK-07C General First-aid Kit is made of aeronautic-aluminum materials as well as broken-proof, shake-proof and antioxidation. With striking Logo which easy to find it's function, safety key, soft gallus and tote for your different way to carriage. It is totally used nontoxic green material inside, while adjustable-clapboard can fix the accessories more fixed. A reasonable layout with transparency medicine-box or organic glass to keep inside instrument laid in order when closing kits. It is very convenient to use in places, especially for long-distance remotion.

Standard Configuration:

No.	Description	Specification	Qty
1	Oxygen Cylinder	1.4 Litre	1set
2	Manual Resuscitator	PVC for Adult	1set
3	Sphygmomanometer	*	1set
4	Stethoscope	*	1set
5	Thermometer	mercury	1
6	Manual Suction	R Type	1set
7	Aerate Equipment for Throat	Size3.4.7.8(disposable)	4
8	Mouth-to-mouth breathing mask	*	2
9	Anaesthetic Laryngoscope	curved(L,M,S)	1set
10	Equipment for opening Mouth	*	1
11	Tongue Forceps	18cm	1
12	Lamina for Pressing Lingua	Metal	1
13	Plexor	*	1
14	Endotracheal Tube	Size8,9(disposable)	4
15	Catheter	Size 14,16	2
16	Gauze Bandage	10*500cm	4
17	Cravat	100*100*140cm	2
18	Compressed Gauze	50*80cm	2
19	Medicinal Gauze Pieces	7.5*7.5cm	10
20	First-aid Carpet	130*210cm	1
21	Rolling Plywood	92*11cm	2
22	Tourniquet	*	1
23	Breathe-freely Adhesive Tape	1.25*200cm	2
24	Medicinal Gloves	Size 7 or 7.5	1
25	Alcohol Cotton	5*5cm	10
26	Iodine Cotton Swabs	5pcs/pac.	4
27	Scissor for Operation	12.5cm	1
28	Dressing Clamp	12.5cm	1
29	Flashlight	*	1
30	First-aid Manual	*	1



EFK-07F Aluminum First-aid Kit for Resuscitation

Size: 44*32*20cm

The EFK-07F Aluminum First-aid Kit for Resuscitation is made of aeronautic-aluminum materials as well as broken-proof, shake-proof and antioxidation. With striking Logo which easy to find it's function, safety key, soft gallus and tote for your different way to carriage. It is totally used nontoxic green material inside, while adjustable-clapboard can fix the accessories more fixed. A reasonable layout with transparency medicine-box or organic glass to keep inside instrument laid in order when closing kits. It is very convenient to use in places, especially for long-distance remotion.

Standard configuration:

No.	Description	Specification	Qty
1	Manual Resuscitator	PVC for Adult	1set
2	Mouth-to-mouth breathing mask*		2
3	Aerate Equipment for Throat	Size3.4.7.8(disposable)	4
4	Sphygmomanometer	*	1set
5	Stethoscope	*	1set
6	Thermometer	mercury	1
7	Manual Suction	R Type	1set
8	Anaesthetic Laryngoscope	curved(L,M,S)	1set
9	Equipment for opening Mouth	*	1
10	Tongue Forceps	*	1
11	Lamina for Pressing Lingua	Metal	1
12	Endotracheal Tube	Size3.4.7.8(disposable)	4
13	Fixer for Endotracheal Tube	*	1
14	Compressed Gauze	50*80cm	2
15	Medicinal Gauze Pieces	7.5*7.5cm	10
16	Breathe-freely Adhesive Tape	1.25*200cm	2
17	Medicinal Gloves	Size 7 or 7.5	1
18	Alcohol Cotton	*	10
19	Iodine Cotton Swabs	5pcs/pac.	4
20	Scissor for Operation	12.5cm	1
21	Dressing Clamp	12.5cm	1
22	Flashlight	*	1
23	Tourniquet	*	1
24	First-aid Manual	*	1



EFK-07W Aluminum First-aid Kit for Trauma

Size: 55.7*39*22cm

The EFK-07W Trauma Aluminum First-aid Kit is made of aeronautic-aluminum materials as well as broken-proof, shake-proof and antioxidation. With striking Logo which easy to find it's function, safety key, soft gallus and tote for your different way to carriage. It is totally used nontoxic green material inside, while adjustable-clapboard can fix the accessories more fixed. A reasonable layout with transparency medicine-box or organic glass to keep inside instrument laid in order when closing kits. It is very convenient to use in places, especially for long-distance remotion.

Standard configuration:

No.	Description	Specification	Qty
1	Manual Resuscitator	PVC for Adult	1set
2	Sphygmomanometer	*	1set
3	Stethoscope	*	1set
4	Thermometer	mercury	1
5	Equipment for opening Mouth	*	1
6	Tongue Forceps	*	1
7	Lamina for Pressing Lingua	Metal	1
8	Flashlight	*	1
9	Scissor for Operation	12.5cm	1
10	Dressing Clamp	12.5cm	1
11	Gauze Bandage	10*500cm	4
12	Cravat	91*91*130cm	2
13	Elasticity Bandage	5*450cm	2
14	Compressed Gauze	50*80cm	2
15	Medicinal Gauze Pieces	7.5*7.5cm	10
16	Rolling Plywood	92*11cm	2
17	Tourniquet	*	1
18	Breathe-freely Adhesive Tape	1.25*200cm	4
19	Medicinal Gloves	Size 7 or 7.5	1
20	Alcohol Cotton	*	10
21	Iodine Cotton Swabs	5pcs/pac.	4
22	First-aid Carpet	130*210cm	1
23	Manual Suction	R Type	1set
24	Neck Holder	*	1
25	Anaesthetic Laryngoscope	curved(L,M,S)	1set
26	Sewing Bundle for Trauma	*	1set
27	First-aid Manual	*	1



EFK-07Z Aluminum First-aid Kit for Resuscitation

Specification: 55.7*39*23.8cm

The EFK-07Z Integrated Aluminum First-aid Kit is made of aeronautic-aluminum materials as well as broken-proof, shake-proof and antioxidation. With striking Logo which easy to find it's function, safety key, soft gallus and tote for your different way to carriage. It is totally used nontoxic green material inside, while adjustable-clapboard can fix the accessories more fixed. A reasonable layout with transparency medicine-box or organic glass to keep inside instrument laid in order when closing kits. It is very convenient to use in places, especially for long-distance remotion.

Standard configuration:

No.	Description	Specification	Qty
1	Oxygen Cylinder	2 Litre	1set
2	Manual Resuscitator	PVC for Adult	1set
3	Sphygmomanometer	*	1set
4	Stethoscope	*	1set
5	Thermometer	mercury	1
6	Manual Suction	R Type	1set
7	Anaesthetic Laryngoscope	curved(L,M,S)	1set
8	Equipment for opening Mouth	*	1
9	Tongue Forceps	*	1
10	Lamina for Pressing Lingua	Metal	1
11	Plexor	*	1
12	Endotracheal Tube	Size3.4.7.8(disposable)	4
13	Gauze Bandage	10*500cm	4
14	Cravat	91*91*130cm	2
15	Compressed Gauze	50*80cm	2
16	Medicinal Gauze Pieces	7.5*7.5cm	10
17	First-aid Carpet	130*210cm	1
18	Rolling Plywood	92*11cm	2
19	Tourniquet	*	1
20	Breathe-freely Adhesive Tape	1.25*200cm	4
21	Medicinal Gloves	Size 7 or 7.5	1
22	Alcohol Cotton	*	10
23	Iodine Cotton Swabs	5pcs/pac.	4
24	Scissor for Operation	12.5cm	1
25	Dressing Clamp	12.5cm	1
26	Flashlight	*	1
27	Sewing Bundle for Trauma	*	1set
28	First-aid Manual	*	1



EFK-07ZT Integrated Aluminum First-aid Kit with Trolley

Specification: 55.7*39*23.8cm

The EFK-07ZT Integrated Aluminum First-aid Kit is made of aeronautic-aluminum materials as well as broken-proof, shake-proof and antioxidation. With striking Logo which easy to find it's function, safety key, soft gallus and tote for your different way to carriage. It is totally used nontoxic green material inside, while adjustable-clapboard can fix the accessories more fixed. A reasonable layout with transparency medicine-box or organic glass to keep inside instrument laid in order when closing kits. It is very convenient to use in places, especially for long-distance remotion.

Standard configuration:

No.	Description	Specification	Qty
1	Oxygen Cylinder	2 Litre	1set
2	Manual Resuscitator	PVC for Adult	1set
3	Sphygmomanometer	*	1set
4	Stethoscope	*	1set
5	Surgical Blade	Size 20	2pcs
6	Manual Suction	R Type	1set
7	Anaesthetic Laryngoscope	curved(L,M,S)	1set
8	Equipment for Opening Mouth	*	1
9	Tongue Forceps	*	1
10	Lamina for Pressing Lingua	Metal	1
11	Plexor	*	1
12	Endotracheal Tube	Size3.4.7.8(disposable)	4
13	Gauze Bandage	10*500cm	4
14	Cravat	91*91*130cm	2
15	Compressed Gauze	50*80cm	2
16	Medicinal Gauze Pieces	7.5*7.5cm	10
17	First-aid Carpet	130*210cm	1
18	Plywood	foldable	1
19	Tourniquet	*	1
20	Breathe-freely Adhesive Tape	1.25*200cm	4
21	Medicinal Gloves	Size 7 or 7.5	1
22	Alcohol Cotton	*	10
23	Iodine Cotton Swabs	5pcs/pac.	4
24	Scissor for Operation	12.5cm	1
25	Dressing Clamp	12.5cm	1
26	Flashlight	*	1
27	Blade Holder	Size 4	1
28	First-aid Manual	*	1



EFK001-CZ1 General First-aid Kit

Size: 32*23*20cm

The EFK001-CZ1 General First-aid Kit is made of materials totally green as well as water-proof and nontoxic. It is portable and light for domiciliary treatment, easy for cleaning and has characters of water-proof, dusk-proof, quake-proof that can be used in places of small-size clinics ,mid-size clinics, families, large and mid-size buses, cars , tourism teams and community hospitals.

Standard Configuration:

No.	Description	Specification	Qty
1	Sphygmomanometer	*	1set
2	Stethoscope	*	1set
3	Lamina for Pressing Lingua	disposable	1
4	Medicinal Scissor	12.5cm	1
5	Flashlight	*	1
6	Dressing Clamp	12.5cm	1
7	Alcohol Cotton	*	10
8	Iodine Cotton Swabs	5pcs/pac.	4
9	Gauze Bandage	10*500cm	6
10	Medicinal Gauze Pieces	7.5*7.5cm	10
11	Cravat	91*91*130cm	2
12	Adhesive Plaster	200*1.25cm	2
13	Compressed Gauze	50*80cm	2
14	Tourniquet	latex	1
15	First-aid Manual	*	1



Standard configuration:

Blood-pressure meter, stethoscope, lamina for pressing lingua, medicine scissors, electric torch, dressing nipper, alcohol cotton, iodine cotton swabs, gauze bandage, medicine gauze pieces, cravat, adhesive plaster, compressed gauze, stanching strip, latex.

Configuration recommended:

Aluminum box, transfusion utensil,injector,10% glucose liquid,0.9% sodium chloride liquid.

EFK001-CZ2 General First-aid Kit

Size: 62*36*30cm

This model is made of materials totally green as well as water-proof and nontoxic. It is easy for cleaning and has characters of water-proof ,dusk-proof and quake-proof. There are diagnose instrument, transfusion instrument, oxygen equipment and resuscitation equipment ,etc. in addition , it is possible to take more than 40 kinds of first-aid syringe and medicine as well as common drugs. It fits well for hospitals , ambulances and large-size clinic and is considered as an ideal medicine bag for doctors who visiting the patients outside.

Standard Configuration:

No.	Description	Specification	Qty
1	Oxygen Cylinder	2 Litre	1set
2	Air Chamber for Resuscitation	PVC for Adult	1set
3	Manual Aspirator	*	1set
4	Sphygmomanometer	*	1set
5	Stethoscope	*	1set
6	Lamina for Pressing Lingua	*	1
7	Equipment for opening Mouth	*	1
8	Tongue Forceps	*	1
9	Thermometer	mercury	1
10	Plexor	*	1
11	Flashlight	*	1
12	Scissor	12.5cm	1
13	Dressing Clamp	12.5cm	1
14	Cervix Fixer	*	1
15	Rolling plywood	92*11cm	2
16	Adhesive Plaster	200*1.25cm	2
17	Gauze Bandage	10*500cm	4
18	Cravat	91*91*130cm	2
19	Compressed Gauze	50*80cm	2
20	Medicinal Gauze Pieces	7.5*7.5cm	10
21	Alcohol Cotton	*	10
22	Iodine Cotton Swabs	5pcs/pac.	4
23	Medicinal Gloves	*	1
24	Tourniquet	latex	1
25	First-aid Manual	*	1



Recommended selection or self-preparation:

Disposable catheter, disposable drainage bag, disposable injector, transparent plaster, trauma dressing, all kinds of syringes and medicine, iodine, styptic paster, compound sodium chloride injection, 0.9 sodium chloride injection, 10% glucose liquid, transfusion equipment, vein acicula, acupuncture acicula.

Common first-aid syringe and medicine and oral administrations listed for recommendations, you can prepare according to each usage. (Adrenalin injection, Isopropyl adrenalin injection, Lobeline injection, injection, Coramine, Dopamine, Reserpine, Furosemide, Cedilan-d Lanatoside

C, Antipsychotic, Atropin, Dicyanine Etamsylate, Dexamethasone, Gentamycin, 654-2 injection, Aninophylline injection, 5% Sodium Bicarbonate, 50% glucose liquid.

Common oral administrations:

Nifedipine, Isosorbide Dinitrate, Glonoin, Antipsychotic.

EFK001-CZ3 General First-aid Kit

Size: 43*26*35cm

The EFK001-CZ3 General First-aid Kit is made up of materials totally green as well as water-proof and nontoxic. It is portable and light for domiciliary treatment, easy for cleaning and has characters of water-proof, dusk-proof that can be used in places of clinics, hospitals, ambulances and community service stations.

Standard configuration:

No.	Description	Specification	Qty
1	Air Chamber for Resuscitation	PVC for Adult	1set
2	Manual Aspitator	*	1set
3	mouth-to-mouth breathing mask	*	1
4	Equipment for opening Mouth	*	1
5	Lamina for Pressing Lingua	disposable	1
6	Tongue Forceps	*	1
7	Flashlight	*	1
8	Medicinal Scissor	12.5cm	1
9	Medicinal Forceps	12.5cm	1
10	Sphygmomanometer	*	1set
11	Stethoscope	*	1set
12	Thermometer	mercury	1
13	Alcohol Cotton	*	10
14	Iodine Cotton Swabs	5pcs/pac.	4
15	Adhesive Plaster	200*1.25cm	2
16	Gauze Bandage	10*500cm	4
17	Cravat	91*91*130cm	2
18	Compressed Gauze	50*80cm	2
19	Medicinal Gauze Pieces	7.5*7.5cm	10
20	Medicinal Glove	*	1
21	First-aid Carpet	*	1
22	Tourniquet	latex	1
23	Rolling plywood	92*11cm	2
24	Medicinal Bottle Opener	Size L & S	2
25	Plexor	*	1
26	First-aid Manual	*	1



Self-preparation recommended:

physiological brine, hydrogen dioxide liquid, asepsis dressing of many kinds of specification, disposable transfusion equipment, injector, ice bag, medicine asepsis Q-tips, disposable catheter, triangular cloth first-aid bag, acupuncture acicula.

EFK001-FS1 First-aid Kit for Resuscitation

Size: 45*22*32cm

This kind of medical bag is made up of materials of water-proof, nontoxic and totally green. It is easy for cleaning and water-proof, dust-proof, quake-proof and has an reasonable layout. It is very convenient to use in places like coastal beauty spot, swimming pools and ferry docks as well as hospitals, first-aid centers, ambulances and community service stations

Standard configuration:

No.	Description	Specification	Qty
1	Air Chamber for Resuscitation	PVC for Adult	1set
2	Manual Aspirator	*	1set
3	Oxygen Cylinder	2 Litre	1set
4	Mouth-to-mouth breathing mask	*	1
5	Sphygmomanometer	*	1set
6	Stethoscope	*	1set
7	Thermometer	mercury	1
8	Laryngoscope	curved(L,M,S)	1set
9	Equipment for opening Mouth	*	1
10	Lamina for Pressing Lingua	disposable	1
11	Flashlight	*	1
12	Tongue Forceps	*	1
13	Medicinal Scissor	12.5cm	1
14	Endotracheal Tube	Size3.4.7.8	4
15	Stanching Nip	12.5cm	1
16	Pharynx Airduct	*	4
17	Medicinal Forceps	12.5cm	1
18	Gauze Bandage	10*500cm	4
19	Cravat	91*91*130cm	2
20	Medicinal Gauze Pieces	7.5*7.5cm	10
21	Compressed Gauze	50*80cm	2
22	Adhesive Plaster	200*1.25cm	2
23	Medicinal Glove	*	1
24	Alcohol Cotton	*	10
25	Iodine Cotton Swabs	5pcs/pac.	4
26	First-aid Manual	*	1



Configuration recommended:

Asepsis wound dressing, wound paster, Q-tips, physiological brine, hydrogen dioxide liquid

EFK001-FS2 First-aid Kit for Resuscitation

Size: 55*32*29cm

This kind of medical bag is made up of materials of water-proof, nontoxic and totally green. It is easy for cleaning and water-proof, dust-proof, quake-proof and has an reasonable layout. It is very convenient to use in places like coastal beauty spot, swimming pools and ferry docks as well as hospitals, first-aid centers, ambulances and community service stations.

Standard configuration:

No.	Description	Specification	Qty
1	Oxygen Cylinder	2 Litre	1set
2	Air Chamber for Resuscitation	PVC for Adult	1set
3	Manual Aspirator	*	1set
4	Endotracheal Tube	Size3.4.7.8	4
5	Sphygmomanometer	*	1set
6	Stethoscope	*	1set
7	Equipment for opening Mouth	*	1
8	Lamina for Pressing Lingua	*	1
9	Tongue Forceps	*	1
10	Flashlight	*	1
11	Medicinal Scissor	12.5cm	1
12	Stanching Nip	12.5cm	1
13	Mouth-to-mouth breathing mask	*	1
14	Medicinal Gauze Pieces	7.5*7.5cm	10
15	Adhesive Plaster	200*1.25cm	2
16	Medicinal Glove	*	1
17	Alcohol Cotton	*	10
18	First-aid Manual	*	1



Configuration recommended:

disposable injector, rubber sucker, prepare laryngoscope if needed.

Attentions:

1. For cleaning :the medicine bag can be cleaned with common oil detergent.
2. For sterilizing:suffocate with 37%-40% formaldehyde liquor.
3. Regularly check the first-aid tools in the medicine bag,and change certain ones if needed.
4. Timely make up for certain instruments and expendables after each usage.
5. The used instruments should be sterilized for the next usage.
6. The configuration above is just for recommendation, they can be configured according to the demand of each usage.
7. For manual absorbing equipment and air chamber for resuscitation, please consult the specification of single product

EFK001-WS1 Trauma First-aid Kit

Size: 44×18×36cm

The medical bag is made of materials of water-proof, nontoxic and totally green. It is easy for cleaning and water-proof, dust-proof, quake-proof and has an reasonable layout. It can be used in places of high-risk. Such as: sports places, building sites and outdoor sports.

Standard configuration:

No.	Description	Specification	Qty
1	Sphygmomanometer	*	1set
2	Stethoscope	*	1set
3	Medicinal Scissor	12.5cm	1
4	Medicinal Forceps	12.5cm	1
5	Thermometer	mercury	1
6	Stanching Nip	12.5cm	1
7	Gauze Bandage	10*500cm	4
8	Cravat	91*91*156cm	2
9	Compressed Gauze	50*80cm	2
10	Medicinal Gauze Pieces	7.5*7.5cm	10
11	Alcohol Cotton	*	10
12	Iodine Cotton Swabs	5pcs/pac.	4
13	First-aid Carpet	130*210cm	1
14	Medicinal Gloves		1
15	Tourniquet	latex	1
16	Cervix Fixer		1
17	Rolling Plywood	92*11cm	2
18	First-aid Manual	*	1



Configuration recommended:

Asepsis dressing of many kinds of specification, wound paster, physiological brine, hydrogen dioxide liquid, cotton, you can prepare some common trauma medicine if needed

EOC-1.4L/2L/3.2L/4.0L/7L/8L/10L/15L/40L Oxygen Cylinder



- For single meter type, It is equipped with pressure meter. Observe the flow humidifier.(1.4L/2L/3.2L/4L)
 - For the double meter type, equipped with pressure meter.and oxygen flowmeter, and manostat device.(1.4L/2L/3.2L/4L)
 - For trolley type, easy to move, long oxygen use time(7L/8L/10L/15L)
 - All size available in steel bottle
 - 2L/10L available in aluminum alloy bottle
- 1.4L -soft bag packing
 2L/3.2L- aluminum alloy case
 4L -leather case/plastic case/aluminum alloy case
 7L/8L/10L/15L-aluminum alloy trolley

Size	Max storage of oxygen	Use time	Referenc e weight	Max oxygen output (single meter)	Max oxygen output (double meter)	Max oxygen pressure
EOC-1.4L	205L	1.5 hour	4kg	≥20L/min	≥4L/min	15MPA
EOC-2L	294L	2 hour	6 kg	≥20L/min	≥4L/min	15MPA
EOC-3.2L	470L	3.5 hour	7 kg	≥20L/min	≥4L/min	15MPA
EOC-4L	588L	4.5 hour	8.5 kg	≥20L/min	≥4L/min	15MPA
EOC-7L	1029L	8 hour	16 kg	≥20L/min	≥4L/min	15MPA
EOC-8L	1176L	9.5 hour	17 kg	≥20L/min	≥4L/min	15MPA
EOC-10L	1470L	11.5 hour	18 kg	≥20L/min	≥4L/min	15MPA
EOC-15L	2205L	26.5 hour	23 kg	≥20L/min	≥4L/min	15MPA
EOC-40L	5880L	47 hour	50 kg			

Note: use time is relative with oxygen pressure and oxygenflow. Above data is estamated accordingly to oxygen 11MPA, oxygenflow 1.5L/min

Soft bag



Leather case



Aluminum case



plastic case



Aluminum alloy trolley

