

AOT400B

Electric Table for Gynaecology and Obstetrics

400B is the most latest electric table that can be used for Obstetrics birthing , Gynecological operation and examination. The tabletop elevation, trendelenburg/reverse trendelenburg, back plate all controlled by both micro touch hand remote control and foot switch. Imported electric motors enable flexible, smooth movements, low noise, high reliability.

All base and column cover is made of 304 medical grade stainless steel. Seamless mattress can be colorfully customized, easy to clean and sterilize.



MEET THE VARIOUS
REQUIREMENTS OF GYNAECOLOGY AND OBSTETRICS



WANROOEMED

WANROOEMED

YOUR IDEAL OPERATING TABLE SUPPLIER
MEET THE VARIOUS REQUIREMENTS OF GYNAECOLOGY AND OBSTETRICS

AOT400B

Electric Table for Gynaecology and Obstetrics



WANROOEMED

YOUR IDEAL OPERATING TABLE SUPPLIER

WANROOE MACHINERY CO., LTD





WANROOEMED YOUR IDEAL OPERATING TABLE SUPPLIER
MEET THE VARIOUS REQUIREMENTS OF GENERAL SURGERY



AOT400B

Electric Table for Gynaecology and Obstetrics

● Examination Position



● Birthing Position



OPTIONAL ACCESSORY


Examination Lamp

Various Colors For Mattress



AOT400B

Electric Table for Gynaecology and Obstetrics
TABLETOP MOVEMENTS



STANDARD ACCESSORIES

Mattress	1 Set
Arm Support	1 Pair
Handle	1 Pair
Fixing Clamp	1 Pair
Knee Crutch	1 Pair
Pedal	1 Pair
Leg Support	1 Piece
Waste Basin	1 Piece
Hand Remote	1 Piece
Foot Switch	1 Piece

MAIN SPECIFICATION

400B Electric Table for Gynaecology and Obstetrics	
Tabletop Length/Width	1780/600mm
Tabletop Up/Down	880/680mm
Tabletop Trendelenburg/Reverse Trendelenburg	8°/20°
Back Plate Up/Down	45°/7°
Voltage/Frequency/Capacity	220V/50Hz/1.0kw
Emergency Battery	Optional

MOVEMENTS DIAGRAM OF AOT400B

