

江苏如翼医疗科技有限公司
JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.
Add: Room1311,1313, Guotai New World Plaza, No19,
East Renmin Road,Zhangjiagang, Jiangsu, China.
Tel:+86 512 58810663
Whatsapp/Mobile:+ 86-18962287398;15962399164
Email:info@wanrooemed.com
Http://wanrooemed.com



www.wanrooemed.com

Electric Wheelchairs



江苏如翼医疗科技有限公司
JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Contents

electric
wheelchairs & scooter

- 01 ● BWHE 301
- 02 ● BWHE301X
- 03 ● BWHE301L
- 04 ● BWHE302
- 05 ● BWHE303
- 06 ● BWHE1801
- 07 ● BWHE1801LX
- 08 ● BWHE1802
- 09 ● BWHE1803
- 10 ● BWHE1002
- 11 ● BWHE1002X
- 12 ● BWHE1003
- 13 ● BWHE501
- 14 ● BWHE501E
- 15 ● BWHE501PE
- 16 ● BWHE502
- 17 ● BWHE503
- 18 ● BWHE601
- 19 ● BWHE602
- 20 ● BWHE603
- 21 ● BWHE901
- 22 ● BWHE2401
- 23 ● BWHEB802
- 24 ● BWHE 803
- 25 ● BWHE201,BWHE02-2201
- 26 ● Certificates



● Foldable Backrest

Push the button, and fold the backrest easily to reduce the space.

● Intelligent control

Ingenious, simple, accurate, full-intelligent, easy to operate.

● Detachable Footrest

Adjustable and detachable. Good quality foot bands, comfortable and safe.

● Flip-up Armrest

Humanization design. Flip-up armrest, easy for users to up and down the wheelchair. Inside is filled with sponge, soft and comfortable, deduct hand pressure effectively.

● Chilwee Battery

Chinese battery industry leading brand.

● Motor

250W DC motor, stable performance.

● Mode switch

Unique safety design. Easy to switch between manual and electric mode.

● Anti-tipper

Avoid wheelchair toppling over. Casters can be easily assembled and disassembled.

electric wheelchairs

BWHE301



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 12Ah/20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: manual brake
- ◉ Frame material: steel
- ◉ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
40/46	1090	685	930	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	420	8	22	100

BWHE301X



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 12Ah/20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
40/46	1090	685	930	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	520	8	22	100

electric wheelchairs

BWHE301L



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 12Ah/20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: manual brake
- ◉ Frame material: aluminium alloy
- ◉ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
36/41	1030	685	955	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	415	8	22	100

BWHE302



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 12Ah/20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: manual brake
- ◉ Frame material: steel
- ◉ Footrest/backrest adjustable angle: 60°/80°
- ◉ Frame color: silver

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
46/52	1110	670	930	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	420	8	22	100

electric wheelchairs

BWHE303



Parameters

- ⊙ Driving range: 15-20km
- ⊙ Max speed: 6km/h
- ⊙ Motor: DC 250W *2Pcs (Jerry)
- ⊙ Battery: 24V 12Ah/20Ah
- ⊙ Control: WANROOE
- ⊙ Parking brake: manual brake
- ⊙ Frame material: steel
- ⊙ Footrest/backrest adjustable angle: 60°/45°
- ⊙ Frame color: silver

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
48/53	1100	685	930	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	420	8	22	100

BWHE1801



Parameters

- ⊙ Driving range: 15-20km
- ⊙ Max speed: 6km/h
- ⊙ Motor: DC 250W *2Pcs (Jerry)
- ⊙ Battery: 24V 12Ah/20Ah
- ⊙ Control: WANROOE
- ⊙ Parking brake: manual brake
- ⊙ Frame material: steel
- ⊙ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
44/49	1110	635	965	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	525	380	10	18	100

electric wheelchairs

BWHE1801LX



Parameters

- ⦿ Driving range: 15-20km
- ⦿ Max speed: 6km/h
- ⦿ Motor: DC 250W *2Pcs (Jerry)
- ⦿ Battery: 24V 12Ah/20Ah
- ⦿ Control: WANROOE
- ⦿ Parking brake: electromagnetic brake
- ⦿ Frame material: aluminium alloy
- ⦿ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
37/42	1030	685	955	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	480	8	16	100

BWHE1802



Parameters

- ⦿ Driving range: 15-20km
- ⦿ Max speed: 6km/h
- ⦿ Motor: DC 250W *2Pcs (Jerry)
- ⦿ Battery: 24V 12Ah/20Ah
- ⦿ Control: WANROOE
- ⦿ Parking brake: manual brake
- ⦿ Frame material: steel
- ⦿ Footrest/backrest adjustable angle: 60°/80°
- ⦿ Frame color: silver

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
49/55	1100	635	980	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
400	530	380	10	18	100

electric wheelchairs

BWHE1803



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 12Ah/20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: manual brake
- ◉ Frame material: steel
- ◉ Footrest/backrest adjustable angle: 60°/45°
- ◉ Frame color: silver

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
51/56	1110	635	980	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	525	380	10	18	100

BWHE1002



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: manual brake
- ◉ Frame material: steel
- ◉ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
46	1140	640	935	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	370	8	12	100

electric wheelchairs

BWHE1002X



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
46	1140	640	935	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	470	8	12	100

BWHE1003



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Jerry)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: aluminium alloy
- ◉ Frame color: gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
28	980	600	950	445	460
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
460	460	N/A	8	12	100













electric wheelchairs

BWHE501



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE , PG, Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Frame color: silver, blue, black, red, gray













 N.W. (kg)	 overall length(mm)	 overall width(mm)	 overall height(mm)	 seat width(mm)	 seat depth(mm)
48	1170	620	940	450	435
 backrest height(mm)	 seat height (mm)	 folded width(mm)	 front wheel(inch)	 rear wheel(inch)	 loading capacity(kg)
395	500	370	8	12	120

BWHE501E



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE , PG, Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Footrest/backrest adjustable angle: 60°/80°
- ◉ Frame color: silver

 N.W. (kg)	 overall length(mm)	 overall width(mm)	 overall height(mm)	 seat width(mm)	 seat depth(mm)
55	1190	620	940	450	435
 backrest height(mm)	 seat height (mm)	 folded width(mm)	 front wheel(inch)	 rear wheel(inch)	 loading capacity(kg)
395	500	370	8	12	120

electric wheelchairs

BWHE501PE



Parameters

- ◎ Driving range: 15-20km
- ◎ Max speed: 6km/h
- ◎ Motor: DC 250W *2Pcs (Taiwan)
- ◎ Battery: 24V 20Ah
- ◎ Control: WANROOE , PG, Dynamic
- ◎ Parking brake: electromagnetic brake
- ◎ Frame material: steel
- ◎ Footrest/backrest adjustable angle: 60°/45°
- ◎ Frame color: silver

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
57	1215	620	940	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	370	8	12	120

BWHE502



Parameters

- ◎ Driving range: 15-20km
- ◎ Max speed: 6km/h
- ◎ Motor: DC 250W *2Pcs (Taiwan)
- ◎ Battery: 24V 20Ah
- ◎ Control: WANROOE , PG, Dynamic
- ◎ Parking brake: electromagnetic brake
- ◎ Frame material: steel
- ◎ Frame color: silver, blue, black, red, gray

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
51	1170	620	940	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	500	370	8	22	120

electric wheelchairs

BWHE503



Parameters

- ⦿ **Driving range:** 15-20km
- ⦿ **Max speed:** 6km/h
- ⦿ **Motor:** DC 250W *2Pcs (Taiwan)
- ⦿ **Battery:** 24V 20Ah
- ⦿ **Control:** WANROOE , PG, Dynamic
- ⦿ **Parking brake:** electromagnetic brake
- ⦿ **Frame material:** aluminium alloy
- ⦿ **Frame color:** silver, gray, blue flame, red flame

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
43	1080	655	945	445	425
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
410	505	410	8	22	120

BWHE601



Parameters

- ⦿ **Driving range:** 15-20km
- ⦿ **Max speed:** 6km/h
- ⦿ **Motor:** DC 250W *2Pcs (Taiwan)
- ⦿ **Battery:** 24V 20Ah
- ⦿ **Control:** WANROOE , PG, Dynamic
- ⦿ **Parking brake:** electromagnetic brake
- ⦿ **Frame material:** steel
- ⦿ **Backrest adjustable angle:** 45°
- ⦿ **Frame color:** red

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
61	1010	600	1075	450	450
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
480	490	N/A	8	12	120

electric wheelchairs

BWHE602



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE , PG, Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: aluminium alloy
- ◉ Frame color: silver, gray, blue flame, red flame

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
43	1060	610	940	450	430
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
395	490	320	8	12	120

BWHE603



Parameters

- ◉ Driving range: 20-30km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 32Ah
- ◉ Control: WANROOE , PG, Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Backrest adjustable angle: 45°
- ◉ Frame color: black

N.W. (kg)	overall length(mm)	overall width(mm)	overall height(mm)	seat width(mm)	seat depth(mm)
74	1010	640	950	450	435
backrest height(mm)	seat height (mm)	folded width(mm)	front wheel(inch)	rear wheel(inch)	loading capacity(kg)
390	425	N/A	8	12	120

electric wheelchairs

BWHE901



Parameters

- ◉ Driving range: 20-30km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 32Ah
- ◉ Control: Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: steel
- ◉ Footrest/backrest adjustable angle: 80°/60°
- ◉ Frame color: black

N.W. (kg)	overall length (mm)	overall width (mm)	overall height (mm)	seat width (mm)	seat depth (mm)
114	1110	630	1280	450	500
backrest height (mm)	seat height (mm)	folded width (mm)	front wheel (inch)	rear wheel (inch)	loading capacity (kg)
650	510	N/A	12	8	120

BWHE2401



Parameters

- ◉ Driving range: 15-20km
- ◉ Max speed: 6km/h
- ◉ Motor: DC 250W *2Pcs (Taiwan)
- ◉ Battery: 24V 20Ah
- ◉ Control: WANROOE , PG, Dynamic
- ◉ Parking brake: electromagnetic brake
- ◉ Frame material: aluminium alloy
- ◉ Frame color: silver, gray, blue flame, red flame

N.W. (kg)	overall length (mm)	overall width (mm)	overall height (mm)	seat width (mm)	seat depth (mm)
44	1080	665	940	445	430
backrest height (mm)	seat height (mm)	folded width (mm)	front wheel (inch)	rear wheel (inch)	loading capacity (kg)
410	490	430	8	24	120

Mobility Scooter

BWHE802



Parameters

- ⦿ Driving range: 15-20km
- ⦿ Max speed: 6km/h
- ⦿ Motor: DC 180W
- ⦿ Battery: 24V 20Ah
- ⦿ Control: Dynamic
- ⦿ Parking brake: electromagnetic brake
- ⦿ Frame color: red

N.W. (kg)	overall length (mm)	overall width (mm)	overall height (mm)	seat width (mm)	seat depth (mm)
55	1180	520	990	420	390
backrest height (mm)	seat height (mm)	folded width (mm)	front wheel (inch)	rear wheel (inch)	loading capacity (kg)
380	565	N/A	9	9	120

BWHE803



Parameters

- ⦿ Driving range: 15-20km
- ⦿ Max speed: 6km/h
- ⦿ Motor: DC 180W
- ⦿ Battery: 24V 20Ah
- ⦿ Control: Dynamic
- ⦿ Parking brake: electromagnetic brake
- ⦿ Frame color: red













N.W. (kg)	overall length (mm)	overall width (mm)	overall height (mm)	seat width (mm)	seat depth (mm)
55	1180	520	990	420	390
backrest height (mm)	seat height (mm)	folded width (mm)	front wheel (inch)	rear wheel (inch)	loading capacity (kg)
380	565	N/A	9	9	120

wheelchairs

BWHE201

- Aluminum frame
- Flip-up armrest
- Quick detachable footrest
- Unique anti-tipper
- Quick release rear spoke wheel with aluminum hand rim















 N.W. (kg)	 overall length(mm)	 overall width(mm)	 overall height(mm)	 seat width(mm)	 seat depth(mm)
16	1030	630	945	450	435
 backrest height(mm)	 seat height (mm)	 folded width(mm)	 front wheel(inch)	 rear wheel(inch)	 loading capacity(kg)
395	500	350	8	24	120

BWHE02-2201

- Aluminum frame
- Low seat height, stable gravity
- Three layers cushion
- User-friendly design, comfortable experience



 N.W. (kg)	 overall length(mm)	 overall width(mm)	 overall height(mm)	 seat width(mm)	 seat depth(mm)
12	1040	580	885	400	410
 backrest height(mm)	 seat height (mm)	 folded width(mm)	 front wheel(inch)	 rear wheel(inch)	 loading capacity(kg)
450	384	310	5	22	120

Certificate

