



江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19,  
East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663

Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanroomed.com

www.wanroomed.com www.rooemed.com

# BATTERY-POWERED ELECTRIC TELEHEALTH CART

The best choice for nursing applications



# BWT-007 SERIES

THE BEST CHOICE FOR NURSING APPLICATIONS!

## INTRODUCTION

BWT-007 series 采用模块化设计以满足文档、药物分配、远程医疗和医疗设备应用。结合了紧凑、重量更轻的设计、优秀的卓越的人体工程学，使其成为护理应用的最佳选择。

非常适合用于检查室、手术室、急诊室、放射科 (PACS)、重症监护 (ICU/C-CU) 单元、PACU 医疗外科区域和放射科 (PACS)。

The BWT-007 series is designed to be modular for documentation, drug distribution, telemedicine and medical device applications. The combination of compact, lightweight design and excellent ergonomics makes it the best choice for nursing applications.

Ideal for examination rooms, operating rooms, emergency rooms, radiology (PACS), intensive care (ICU/CCU) units, PACU Medical surgical areas, and radiology (PACS).

# BWT-007

## PC Computer System

Equipped with windows-10, function customization.

## Platform

Provide enough space for the files, scanners, and other objects you need.

## Base

High strength frame base, more stable chassis.

## Power Supply

Display real-time power.

## Power Module

DC output, reliable, safe.

## 4-inch Castors

You can move fully between rooms, or lock casters to secure.



### Lift Button

With motorized lift system to provide an easy and quick height adjustment.



### Monitor Bracket

Height adjustable, and available for tilt position and horizontal rotation. VESA hole: 75\*75mm, 100\*100 mm



### Keyboard Tray

Pull-out design stores keyboard and Bluetooth mouse.



### Lift Column

Aluminum alloy profile lifting column, beautiful and durable.



## Customizable Storage Drawer

Provide extra space for your important documents and equipment.

# BWT-007 Plus

## Camera

USB 3.0 HD conference camera.

## PC Computer System

Equipped with windows-10, function customization.

## Platform

Provide enough space for the files, scanners, and other objects you need.

## Customizable Storage Drawer

Provide extra space for your important documents and equipment.

## Optional Loudspeaker

Real sound, noise suppression, excellent performance.

## Power Supply

Display real-time power.

## Base

High strength frame base, more stable chassis.

## 4-inch Castors

You can move fully between rooms, or lock casters to secure.

## Power Module

DC output, reliable, safe.

## TECHNICAL PARAMETERS

Model	BWT-007 Series	
Overall Size	600*475*1166~1516mm	
Platform	Size	500mm*475mm
	Adjustable Height	795mm-1145mm
Bracket	Monitor Size	≤22 inch
	Loading	6.6-26.4lbs (3-12 kg)
	Lift range	90mm
	Adjustable tilt angle	0° to 15°
	Horizontal rotation angle	-90° to 90°
	Drawer	Small Drawer
Medium Drawer		390*310*160mm
Base	Size	484*408mm

## HIGHLIGHTS

1. Pearl white paint base, easy to clean design helps prevent infection.
2. The height-adjustable workspace makes users comfortable while sitting or standing.
3. Provide battery compartment and clear power display.
4. Organized cable management provides a neat look.
5. Helps improve workflow efficiency and upgrade patient care.
6. DC power system promotes universal compatibility and long battery run time.
7. Diversified video system to realize real-time interactive audio, video communication and sharing of medical resources.
8. The audio system can make the medical staff communicate more smoothly, and at the same time meet the needs of certain professional fields.



### USB 3.0 HD conference camera

Size	220×173×117mm
USB3.0 Video interface	USB3.0
Horizontal rotation angle	-170° ~ +170°
Elevation angle	-35° ~ +90°
Input voltage	DC12V
Image sensor	1/2.8 inch high quality, 5 million HD CMOS sensor
Resolution ratio	1920x1080p/30/25, 1280x720p60/50/30/25
Multiplicity	12X optical times + 12X electron times
Focal length	f=3.92mm (near)~47.32mm (far)
Horizontal viewing angle	6.3° (near) ~ 72.5° (far)
Camera Tray	190×124mm

## ACCESSORIES



### Computer

PC Computer System	CPU	Intelcore i5-10200H/8 G DDR3/256G SSD
	System	Windows 10 & Windows 7
Display	Display size	21.5-inch TFT LCD with 1920x1080 pixels
	Optional	Dual Monitor

## POWER SUPPLY(DC OUTPUT)



### Battery

Input	AC 100V-240V 60Hz/50Hz
Output	DC-12V 6A DC-24V 6A
Output Power	200W
Weight	5.5kg
Optional	UPS



Optional: UPS power system AC output

**OPTIONAL**



Storage Basket



2L Needle Disposal Holder



Multilayer Drawer



Scanner



Scanner Holder




Printer Tray



Mouse Tray



Dual Screen Bracket



Workstation Cart With Multi-drawer

**Specifications**

Accessories	Size
Small Drawer	390*310*55mm
Medium Drawer	390*310*160mm
Printer Tray	346*229mm
Mouse Tray	185*180mm
Keyboard Drawer	270*400mm
Scanner Holder	120*100mm
2L Needle Disposal Holder	Φ140*100mm
Storage Basket	200*128*100mm



## Integrated Mobile Medical Carts Solutions

■ BWT-007 Mobile Powered Workstation Cart

■ BWT-007 Plus Mobile Powered Workstation Cart With Camera



One Drawer



Two Drawers



Two Drawers



Dual Screens