


ROOE 如翼

STAIR CHAIR


爬楼椅



we offer stair chair solutions
that enable safe
evacuation and climbing

我们提供安全疏散以及
爬楼的楼梯椅解决方案

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.
江苏如翼医疗科技有限公司



江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19,
East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663

Whatsapp/Mobile:+ 86-18962287398;15962399164


Email:info@wanrooemed.com

www.wanrooems.com www.wanrooemed.com

INTRODUCTION

我们精心设计的爬楼椅为行动不便的人爬楼带来较大的便利，它的高度的通用性，可以很容易地作为消防疏散椅、紧急疏散椅，或用于简单的家庭使用。

Our well-designed climbing chair brings great convenience for people with mobility difficulties to climb the stairs. With a high degree of versatility, it can be easily used as a fire evacuation chair, emergency evacuation chair, or for simple family use.



we offer stair chair solutions that enable safe evacuation and climbing.

我们提供安全疏散以及爬楼的楼梯椅解决方案 ✓

APPLICATION



Public Building



Office



Retail Establishment



Hospital



Shopping Center



Emergency Service



Health Care Center



Personal Use

Evacuation chairs can be used on almost all types of stairs where upper and lower passages are provided for the safe evacuation of non-moving people.

Our crawler stair climbers are suitable for one-person operation. Made of high-strength aluminum alloy material, Lightweight yet incredibly strong, the chairs can transport persons of up to 180Kg. With flat seat, armrests and padded headrest to provide a comfortable experience.

Castors with brakes for mobility.

All-important rubber tracks help caregivers function on any surface and ensure safe descents. The backrest is designed with two folding handles, allowing the operator multiple grip options for safe transport. Once on the ground, it can be used like a regular wheelchair. Folds into a compact design, making it easier to store and transport. Suitable for indoor and outdoor use, it is economical and requires minimal maintenance.



OPTIONAL

... Colour



Red



Green



Yellow



Blue



White

... Seat



ABS Board



PVC Leather With Foam



PVC Leather



Aluminium Plate

... Other



Dust Cover



Sponge Pad



Wall Bracket



Product Accessories



Curved Handrail



Big Wheels



Fifth wheel



Front Carrying Handles

选配

EA-6FPN

Motorized Stair Climbing Chair(Adjustable Speed)



Casters With Brakes



Front Carrying Handles



Handrail



Folding Structure



Backrest




Seat Belt




Speed Adjustable



Track Switch



Rubber Belt Track



Motorized Stair Climbing Chair

EA-6FPN MINI

Economy Type(Steel Structure,Fixed Speed)



EA-6FPN BASIC

Economy Type(Fixed Speed)



Control Panel



Backrest Adjustment Switch



Seat Folding Switch



Track Folding Structure



Seat Armrest



Economy Type



Castor



Optional Fifth Wheel
(Suitable for narrow stairs)



Cross Grain Track



Track Switch

ROOE 如翼
EA-6FPN PLUS
 Motorized Stair Climbing Chair



Jumbo Honeycomb Wheels With Brakes




Folding Structure

TECHNICAL PARAMETERS

Model	Product Size	Folded Size	Package	N.W/G.W.
EA-6FPN	108*52*120cm	108*52*25cm	116*58*32cm	33kg/46kg
EA-6FPN Basic	108*52*120cm	108*52*25cm	123*58*32cm	33kg/46kg
EA-6FPN Mini	102*57*120cm	110*57*30cm	115*65*37cm	38kg/51kg
EA-6FPN Plus	110*62*120cm	114*62*35cm	115*69*45cm	40kg/55kg



ROOE 如翼
EA-7FPN
Cargo Hand Truck



PRODUCT INTRODUCTION

Electric stair climbing trunks can help operators safely and quickly lift heavy objects up and down stairs. Greatly improve the efficiency of climbing handling, reduce labor costs.

TECHNICAL PARAMETERS

Model	Folded Size	Product Size	Packing Size (1PC)	Loading	N.W	G.W
EA-7FPN	115*50*27cm	80~110*50*110~160cm	116*55*30cm	≤160kg	27kg	40kg

ROOE 如翼
EA-8FPN
 Stair Climbing Wheelchair




PRODUCT INTRODUCTION

An electric stair climber and wheelchair in one product. It can be used directly as an electric wheelchair on flat ground for daily mobility and can be switched to electric mode for stair climbing when you encounter stairs, solving the problem of not being able to go up and down stairs.

FEATURES

1. With two driving modes, easy for climbing and movement on ground.
2. Anti-tipping castor to ensure safety.
3. 4 Key Joystick Controller.
4. Bigger rear tires, stronger stabilization.




TECHNICAL PARAMETERS

Model	Motor	Battery	Folded Size	Product Size	Packing Size(1PC)	Loading	N.W	G.W
EA-8FPN	200W 250W*2	24V/18A	108*64*30cm	108*64*103cm	116*70*36cm	≤181kg	45kg	60kg

EA-6FPNW

Wheelchair Assisted Climbing Chair



PRODUCT INTRODUCTION

This battery operated stair trolley converts any wheelchair into a stair climber.

FEATURES

1. Multi-gear adjustment of aluminum alloy buckle.
2. Wheelchair slot, effectively fix the wheelchair.
3. 82cm length track, can take 3 steps, no bumps when going up and down stairs.
4. Powerful chassis for wheelchair docking.



TECHNICAL PARAMETERS

Model	Motor	Battery	Width Range	Product Size	Packing Size(1PC)	Loading	N.W	G.W
EA-6FPNW	200W	36V/15.6A	40-80cm	108*55*78cm	117*51.5*39cm	≤200kg	36kg	48kg



EA-6FPNW PLUS
Separate Type Climbing Wheelchair



PRODUCT INTRODUCTION

One chair for two uses, can be quickly dismantled. Can be used for electric stair climbing or as a care wheelchair, suitable for use in homes, hotels, etc. for going up and down stairs.

FEATURES

1. Together it is a stair climbing chair, when separated it is a wheelchair.
2. Rear-mounted stair climber easily copes with many environments.
3. Care wheelchair to achieve easy travel.
4. Adjustable thickened headrest, adjustable according to height.

TECHNICAL PARAMETERS

Model	Motor	Battery	Product overall Size	Stair lift size	N.W	Wheelchair size	Wheel chair folded size	N.W
EA-6FPNW 200W Plus	24V/13AH		93~108*56*97~153cm	23~38*40*97~153cm	18.7kg	56*100*99cm	70*53*50cm	17kg



EA-5FPN

Powered Electric Climbing Chair



FEATURES

1. The chair is balanced and does not tip back.
2. Carbon steel material, track undermounted, smooth climbing.
3. One click to place/retract track without any manpower.
4. Armrest can be lifted, easy to fold.
5. Shock absorbing super honeycomb wheel.

PRODUCT INTRODUCTION

Thanks to the uniquely functional climbing track, it is possible to grip the steps for a smooth up and down experience. It is very handy as a stair climbing chair with a powerful motor and high-capacity lithium battery, and as a wheelchair with great mobility independence, we make your life easier.

EA-5FPN PLUS

Powered Electric Climbing Chair With Joystick



FEATURES

1. Two driving modes: wheelchair mode and climbing mode
2. Intelligent controller controls the chair running on the ground and supports a variety of road conditions.

TECHNICAL PARAMETERS

Model	Motor	Battery	Folded Size	Product Size	Packing Size(1PC)	Loading	N.W	G.W
EA-5FPN Plus	200W 250W*2	24V/18A	88*63*48.5cm	96*63*115cm	95*72*54cm	≤181kg	45kg	63kg
EA-5FPN	200W	24V/18A	88*63*48.5cm	96*63*115cm	95*72*54cm	≤181kg	42kg	60kg

ROOE 如翼
EA-6G
Evacuation Chair



PRODUCT INTRODUCTION

Evacuation Chair provides a stair mobility solution when you need it most, reduces health and safety concerns during use and allows for a smooth decent of stairs.

FEATURES

1. Track transmission structure is adopted to facilitate patient decent of stairs.
2. Adjustable head pad and comprehensive patient restraint.
3. Special wheel design provides excellent tracking and stability in confined Spaces and narrow stairwells.
4. Curved handrails reduce waist pressure.
5. The patient can normally be lifted by one person. This reduces health and safety concerns during use and allows for a smooth decent of stairs.
6. It can fold up and be stored away neatly against a wall or in a box. A simple, lightweight and easy to use product.
7. optional: dust cover, wall mounts, front carry handles.



TECHNICAL PARAMETERS

Model	Folded Size	Product Size	Packing Size (1PC)	Loading	N.W	G.W
EA-6G	110*53*20cm	92*53*140cm	116*56*22cm	≤160kg	11kg	13KG

ROOE 如翼
EA-6F
Evacuation Chair



PRODUCT INTRODUCTION

The EA-6F evacuation chair is the ideal equipment for pass through stairs or uneven surfaces with an injured person in case of emergency. The chair has a special track structure which makes it easy to use and safer to move.

Our evacuation chair with 4 wheels, and it enables the patients to travel over most surfaces with minimal effort.

TECHNICAL PARAMETERS

Model	Folded Size	Product Size	Packing Size (1PC)	Loading	N.W	G.W
EA-6F	98*54*35cm	98*54*80cm	110*52*27cm	≤160kg	15kg	17kg

