

江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plazza, No19, East Renmin Road, Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663

Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com

www.wanrooemed.com www.rooemed.com







BLOOD DONOR CHAIR

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.



Electric Two function Blood Donor Chair

BXD100B

The Blood donation chair offers better benefits to individuals such as volunteers who donate blood because it can also be used as a dialysis chair, a medical and rehabilitation chair, and even a bed.

Electric Blood donation chair offers adjustable and flexible sitting surfaces in addition to its already impressive technical features and specifications. Some of the features of Blood Chair include a high quality steel frame, multiple functional options, and maximum load bearing capacity totaling 200Kg.



- 1. Designed according to the principle of human body engineering, luxurious, comfortable and durable. You won't get tired if you sit too long.
- 2. The blood collection chair consists of four main parts: backrest, seat, leg plate and foot pedal.
- 3. It carries out the linkage movement of back, legs and back legs through the electric remote control.
- 4. Wide and comfortable handrails support blood collection.
- 5. Extendable foot pedal boards provide patient comfort and comfort.
- 6. Q235 steel pipe is durable, cushion is made of imported PVC high-quality artificial leather, polyurethane one-time molding.
- 7. Four-inch casters, each wheel can be individually locked and the seat can be fixed in the desired position.
- 8. Standard saltwater rack.
- 9. The chair can be easily used for blood collection and dialysis procedures. 10.Optional: Dining table; Manual CPR.








Seat position

Semi-fowler position

Horizontal position







CPR trendelenburg position

Manual CPR (optional)

Dining table (optional)

roduct model	BXD100B
Chair open size	1900*850*560mm
hair body size	920*850*560mm
Package size	1620*820*600mm/72kg
Seat width	600mm
Arm rest adjustment height	none
	Chair Adjustment Range
Height Adjustable	none
Back adjustment Angle range	-10° ~ 75°±5°
eg plate Adjustment Angle range	-80° ~ 12°±5°
rendelenburg position	-12°
V pole (from clamp to top)	860~1500mm
/ pole (from floor to top)	1260~1900mm
Paper holder and oxygen holder diameter	none
	Material
rame Material	Steel
latress Material	PVC
Vheel	4 Silent Wheels With Brakes
Max loading	200kg
N.W/G.W.	62kg/75kgs
Power Supply	AC220V 50Hz/60Hz



Manual Reclining Blood Draw Chair

BXD100B Basic

Our Reclining Phlebotomy chair can be adjusted by a handle next to the seat for examination of the patient. The blood collection chair uses a gas spring mechanism instead of electric adjustment. Adjustable Phlebotomy Chair's wheels move the chair from one place to another. The maximum safe bearing capacity of 200kg also provides high stability in the horizontal position. Phlebotomy draw chairs are widely used in blood donation centers and medical centers.



- 1. Low carbon steel solid structure, appearance of high strength fiber material and multi-layer baking paint, with artistic sense, modern sense, durability and easy to clean.
- 2. The mattress is made of high-quality medical polyurethane foam forming, which is easy to clean and maintain. Its antistatic and antifouling performance is superior to other similar products.
- 3. The combined backrest and legrest movement promotes a relaxed reclining/lying flat position through gas spring assisted movement.
- 4. Release the handle to quickly adjust the seat position. The release handle can be operated from either side of the recliner.
- 8. Super wide upholstered handrails maximize comfort. Optional dining table board.
- 9. Standard IV pole
- 10.It is an ideal blood collection choice for blood stations, centers and blood donation rooms.
- 11. Human-oriented original products can give blood donors a sense of security and comfort, relieve psychological pressure, and make the whole process relaxed and happy.







Seat position

Semi-fowler position

Backrest adjustment, backleg linkage







Horizontal position

Cpr trendelenburg position

Manual adjustable body position

Product model	BXD100B Basic	
Chair open size	1900*850*560mm	
Chair body size	920*850*560mm	
Package size	1620*820*600mm/72kg	
Seat width	600mm	
Arm rest adjustment height	none	
	Chair Adjustment Range	
Height Adjustable	none	
Back adjustment Angle range	-10° ~ 75°±5°	
Leg plate Adjustment Angle range	-80° ~ 12°±5°	
Trendelenburg position	-12°	
IV pole (from clamp to top)	860~1500mm	
IV pole (from floor to top)	1260~1900mm	
Paper holder and oxygen holder diameter	none	
	Material	
Frame Material	Steel	
Matress Material	PVC	
Wheel	4 Silent Wheels With Brakes	
Max loading	200kg	
N.W/G.W.	62kg/75kgs	



Electric Three function Blood Donor Chair

BXD100B Plus

The adjustable phlebotomy chair has an adjustable seat height and is designed to meet the various healthcare needs of patients for examination, diagnosis and treatment. Electric phlebotomy chairs have a high-density foam cushion that is flexible enough not to cause postural strain when sitting or lying for long periods of time, and PVC leather that is durable and luxuriously comfortable. The dialysis recliner chairs' retractable, foldable footrests provide comfort for patients during dialysis.



- 1. Steel frame coated with epoxy powder with high mechanical strength, chemical and corrosion, beautiful, modern, durable and easy to clean.
- 2. The mattress is made of high-quality medical polyurethane foam forming, which is easy to clean and maintain. Its antistatic and antifouling performance is superior to other similar products.
- 3. With its easy-to-use controller, the seat can stop in any intermediate position and the patient can be positioned comfortably.
- 4. The hand controller adopts easy-to-read button design, can automatically adjust the sitting position or flat bed, and can individually adjust the seat Angle, seat height, backrest Angle, pedal Angle, CPR emergency headdown positioning (or trendlenburg position).
- 5. Four 4-inch silent medical casters with independent brakes.
- 6. Removable headrest The pillow can be adjusted according to the height of each patient.
- 7. Cushioned armrest for comfortable sitting position.
- 8. Electric power control and adjustment of sitting position and foot position to create a safe and comfortable environment for blood donors.
- 9. There are holes for brine rack around, which is convenient to plug and install the brine rack in the plug pipe.
- 10. Ideal for dialysis, blood collection, respiratory and cardiac care.
- 11.Optional: Dining table ; Manual CPR.








Seat position

Horizontal position

CPR Trendelenburg position







Hi-Lo height adjustability

Manual CPR (optional)

Dining table (optional)

Product model	BXD100B Plus	
Chair open size	1900*850*460*830mm	
Chair body size	920*850*460*830mm	
Package size	1620*820*600mm/72kg	
Seat width	600mm	
Arm rest adjustment height	none	
	Chair Adjustment Range	
Height Adjustable	460~830mm	
Back adjustment Angle range	-10° ~ 75°±5°	
Leg plate Adjustment Angle range	-80° ~ 12°±5°	
Trendelenburg position	-12°	
IV pole (from clamp to top)	860~1500mm	
IV pole (from floor to top)	1260~1900mm	
Paper holder and oxygen holder diameter	none	
	Material	
Frame Material	Steel	
Matress Material	PVC	
Wheel	4 Silent Wheels With Brakes	
Max loading	200kg	
N.W/G.W.	62kg/75kgs	
Power Supply	AC220V 50Hz/60Hz	



Electric Three function Blood Donor Chair

BXD100B Pro

The nursing recliner makes patient care, treatment and recovery more comfortable and provides greater convenience for caregivers in a variety of clinical Settings. Helps your patients relax while improving efficiency.



- 1. Fully electric and accurately adjustable blood drawing chair, which can easily locate blood donors.
- 2. Electric remote control adjustment: overall lifting, back lift, leg lift, back leg linkage.
- 3. Double armrest design, the upper armrest can be adjusted left and right, to avoid the blood donor arm sliding and unstable placement.
- 4. Velcro headrest for head and neck comfort during blood donation.
- 5.4 Casters with brakes for easy movement.
- 6. Foldable foot pedals to provide adequate foot support and prevent the patient from sliding down.
- 7. Manual CPR to assist in emergency situations.








Seat position

Backrest adjustment, backleg linkage

Horizontal position







CPR Trendelenburg position

Dining table (optional)

Movable upper handrail

The reclining phlebotomy chai chair can be powered upright, partially reclined or fully reclined. With synchronized headrest and leg support for easy adjustment and quick recovery of positioning. With an upholstered seat and a back with adjustable armrests, patients are sure to stay comfortable. The electric phlebotomy chair keeps patients comfortable and can better accommodate patients who may faint.

Product model	BXD100B PRO
Total Length	1950mm+20mm
Seat width	600mm (total width 900mm including armrest)+20mm
Seat height	580-780mm+20mm
Backrest length	860mm+20mm
Height of armrest	230~290mm+20mm
Castor	4xφ 100mm swivel castors with seperate brakes
Pillow size	410mm*230mm*80mm+20mm
Backrest Adjustment	(-12°~75°)+5°
Largest Adjustment	(-75°~12°)+5°
Weight	70KGS+3KGS
Safe Maximum Load	240KGS