





USKX5000

FULL DIGITAL VETERINARY ULTRASOUND SCANNER



WANROOE
MEDICAL



● Parameter:

- Net weight: 2.2kg.
- 10.4" High resolution color LED monitor, high brightness, high contrast, wide visual, LED back light, image clear and exquisite.
- With 2pcs Large capacity lithium battery(2700MAh), support 200 minutes working time.

● Feature:

- Laptop, slim design, smart and light, convenient for carrying, easier for examining patients at home.
- Menu operation, support 5 customized languages: English, Spanish, Russian, French, Portuguese.
- 10 kinds of Magnification, maximum display depth is 236 mm.
- 100 frames permanent storage, 256 frames cine loop.
- OB measurement: EDD and GA for Bovine, equine, ovine, canine, feline, goat, swine and llama.
- General measurement, Cardiac measurements, Back fat and Lean percentage of swine.
- OB report function.
- USB port, support large capacity USB disk, and color laser printer directly print various image and reports.
- Dual-mode TV output: P A L / N T S C, connecting to thermal printer and ultrasonic working station.



6.5MHz multi-frequency linear rectal probe



6.5MHz multi-frequency micro-convex probe




3.5MHz multi-frequency convex probe



7.5MHz multi-frequency linear probe

● Standard configuration:

Main unit: 1 pc
 6.5MHz linear rectal probe and probe support : 1 pc for each
 Power adapter: 1 pc
 Internal Li-battery : 2 pcs
 Metal case: 1pc
 Battery Charger: 1pc



万容机械有限公司 WANROOE MACHINERY CO., LTD.
 江苏省张家港市人民东路19号国泰新天地广场1311,1313室
 Room1311,1313, Guotai New World Plaza, No19, East Renmin Road,Zhangjiagang,
 Jiangsu, China.
 Tel:+86 512 58810663 Fax:+ 86 512 58810663
 Email:doris@wanrooemed.com