



USKX5600

Full Digital B Mode Ultrasound Scanner



Applicable scope



Applicable Field

USKX5600 is a light weight portable ultrasound well-designed by the export team in Research Center for Ultrasound Medical Applicable Engineering Technology of Jingsu Province, possessing high reliability and performance . It's the ideal option for the pet hospitals, breeding swine farm, animal husbandary and veterinary station, livestock improvement station, stud-farm,dairy farm,marine museum and scientific research institation etc..



USKX5600

Full Digital B Mode Ultrasound Scanner



Reasonable Humanization Design

- 10.4' color LED full screen display, angle adjustment
- Rapidly one-key operation includes one-key storage one-key browse and one-key print
- Large-capacity lithium battery can drive the main unit to work more than 200 minutes
- Internal 4G C disk storage capacity, permanent 5000 frames image, without loss when power off

Powerful Data Management

- Rapidly retrieve the image data management function
- Internal work station multi-plan template, realize process multi medical histories function, quickly finish data, report etc. Manage function
- I-Station integral work station, connect to network, realize one-key print

Abundant Software Package

- OB measurement : GA and EDD for swine, bovine, canine, equine, feline, llama and ovine
- Other measurement: Calculate lean percentage of swine
- OB report function



USKX5600 Full Digital B Mode Ultrasound Scanner

High quality image



Gravidity of dog which is 30 days



Dog's kidney



Dog's bladder

Standard configuration

Main unit

6.5 MHz micro-convex small animal probe

Net cable

Power adapter

Optional configuration

3.5 MHz linear lion muscle probe

3.5 MHz convex probe

6.5 MHz linear rectal probe

7.5 MHz linear probe

Lithium battery 4400 mAh

Charger

Trolley

Vedio recorder (up-897MD)

Enthernet switch

Standoff for back-fat probe

Black handbag



万容机械有限公司

WANROOE MACHINERY CO., LTD.

江苏省张家港市人民东路19号国泰新天地广场1311,1313室

Room1311,1313, Guotai New World Plaza, No19,
East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663

Fax:+ 86 512 58810663

Email:doris@wanrooemed.com



Product depends on goods, this datum is
subject to change without further notice