


BAE507 Five Functions Electric Hospital Bed



Feature:

- Made of superior cold-rolled steel
- One step forming sleeping bed
- Head and foot board, are all made of PP.
- Remote control panel is equipped on the side rails which makes patients operate much easier.
- 4 section bed surface
- 6-bars-Aluminum alloy siderails
- Economic medical castor in 5 inch diameter, two with brakes
- remote controller on the siderails



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com

BAE507 Five Functions Electric Hospital Bed



Standard Specifications

Model number	BAE507
Measurement	Overall size(L*W*H): 2150*950*400-750mm(with 5 inch diameter castor) Patient Touching surface: 1950*830mm
Weight Load	Fully tested sturdy construction capable of taking a maximum user weight of up to 250kg
Function	Backrest max upward angle: 75°
	Footrest max upward angle: 40°
	Height: 400-750mm (castor diameter: 5 inch)
	Trend/Reverse Trend.±12°
Power	220VAC or 110VAC Frequency: 50-60HZ Power Plug: EU type standard if other type, confirm each order

Standard Accessories List

Headboard/Footboard	1 set
Siderails	4 pcs
I.V pole	1 piece
Castor	4 pieces/5 inch
Drainage Hook	2 pieces
Motor	4 pcs

Package: 2080*990*330mm/Carton, Gross weight: 107kgs

20GP: 42 pcs 40GP: 84pcs 40HQ: 96pcs



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com