


# DM-6C Veterinary Anesthesia System

*Friendly Powerful Reliable*

## Application

Designed for safe, simple and controllable inhalation anesthesia for animal.

The Veterinary Anesthesia machine provides rapid control of anesthetic depth with little resistance to respiration from high to low acuity, simple to complex cases, small to big animals.



Flowmeter



Pressure gauge



Breathing bag



Basket



Veterinary Mask  
(Optional)



Cherish your life, Cherish your health!

**WANROOEMED**



# DM-6C Veterinary Anesthesia System

## Trust Points

- Simplicity: Four static casters with self-locking function.
- Option: Wide range application for animals within 200Kgs.
- International standard and advanced technology suitable for wide range use.
- Pressure gauge and flow meter ensure accurate adjusting of parameters.
- Over 2,000 units are installed in more than 200 countries.
- Designed and manufactured by Perlong Medical with over 23 years of experience in veterinary field.

## Features

- High precision flow meter, instantly know the fresh gas flow to the animals.
- Integration breathing circuit design, ensure easy operation and keep tidy.
- Vaporizer with temperature, flow compensation and self-lock function, keep safety anytime.

## Vaporizer



Halothane: 0~5%    Enflurane: 0~5%    Isoflurane: 0~5%    Sevoflurane: 0~8%

1 set vaporizer for standard, can be installed 2 vaporizers in all.

## Specifications

Flow meter	O <sub>2</sub> (0~5L/min)
Rapid oxygen supply	35L/min~75L/min
<b>Monitoring parameters</b>	
Airway pressure	0 cmH <sub>2</sub> O~100cmH <sub>2</sub> O
<b>Working conditions</b>	
Gas source	O <sub>2</sub>
Pressure	280kPa~600kPa

<b>Main component wooden case packing size</b> L 645*W 590*H 630mm	
G.W.	45KG
N.W.	20KG
<b>Accessory wooden case packing size</b> L 650*W 690*H 520mm	
G.W.	40KG
N.W.	20KG

**WANROOEMED**

WWW.ROOEMED.COM

