

NLF-200C CPAP System

Friendly Powerful Reliable

Application

The CPAP System makes a good performance in Neonatal Intensive Care Units (NICU) and other departments. Rooe Medical CPAP System offers you the choices you need in ventilation, monitoring and technique. The expertise in ventilation is based on more than 23 years rich history of developing CPAP solutions that meet your requirements. Great performance is based on high security, high accuracy, high stability and accurate monitoring.




Cherish your life, Cherish your health!



NLF-200C CPAP System

Trust Points

- Designed and manufactured by Rooe Medical with over 23 years experience in this field.
- International standard and advanced technology is suitable for all baby patients around the world.
- 3.5" TFT LCD screen displays the ventilation parameters and alarm information.
- High precision oxygen concentration shows in real time.
- Built-in air and oxygen blender, ensure stable oxygen concentration.
- Advanced power management control technology.



3.5" TFT LCD screen



Shuttle button



Humidifier

Technical Specifications

| | |
|----------------------|---|
| Oxygen concentration | Adjusting range: 21% - 100% Monitoring range: 15% - 100% |
| CPAP | Adjusting range: 0cmH ₂ O - 15cmH ₂ O Monitoring range: 0cmH ₂ O - 60cmH ₂ O |
| Continuous flow | Adjusting range: 2L/min - 20L/min Monitoring range: 2L/min - 20L/min |
| Gas source | O ₂ , AIR |
| Pressure | 280kPa - 600kPa |

| | |
|--------------------------|---|
| Voltage | 100 - 240V |
| Power frequency | 50Hz ± 1Hz |
| Input power | Main component: ≤ 25W Air compressor: ≤ 675W |
| Wooden case packing size | L 670*W 700*H 1410mm |
| G.W. | 90KG |
| N.W. | 57.2KG |

江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19,East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663 Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com Http://wanrooemed.com