

DF-70 PT Training Bed



It is used for fixing patients' different parts so as to prevent following moving when traction treatment.

Size: 196*74*40cm

Weight: 45kg

Load: 135kg

DF-71 PT Training Bed



Used for PT training patient to do exercise on the bed.

Size: 191*125*48cm

Weight: 53kg

Load: 135kg

DF-72 Electric PT Training Bed

Used for PT training patient to do exercise on the bed.

Power: AC220, 50HZ

Size: 199*120*50-80cm

Weight: 127kg

Load: 135kg

DF-73 Electric PT Training Bed

Used for PT training patient to do exercise on the bed.

Power: AC220, 50HZ

Size: 196*123*57-88cm Weight: 134kg Load: 135kg