

DG100 Hydraulic Patient Hoist



DG100 Hydraulic Patient Hoist is suitable for stay in bed for a long term, Patients with inconvenient action or can not move after operation. They use it accompanied by nursing staff, Hospital, nursing center or family can be used.

The main frame is made of high strength steel, the powder of surface plasticsprayingtreatment, The front wheel of the castor is 3 inch universal wheel,the back wheel is 4 inch brake wheel, foot stool can adjust leg width, Sling suspension is safe and reliable, The hydraulic lifting device, the scope is broader.



Specifications

Weight Capacity	150kgs standard
Overall Length	1100mm
Lifting range	740-1775mm
Legs open width	800mm
Ground clearance	120mm
Front twin casters	3"
Rear break casters	4" with brakes
Power	220V,50HZ standard
Pull handle height	1130mm
Charger	Output: 24V DC
Warranty	Limited One (1) Year Warranty
Package	1160*660*250mm, 144pcs/20GP, 324pcs/40GP, 360pcs/40HQ, gross weight: 41kgs, net weight: 38kgs



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com

DG203

Battery-Powered Assist Lift



General Description :

WANROOEMED Battery-Powered Assist Lift is one ideal transport device to provide quick, safe assistance to those who have trouble standing. With multi-safety design, DG203 Battery-Powered Stand-up Lift ensures caregiver and resident security

Features :

- 1.Heavy-gauge steel construction
- 2.Leg open function to mostly close the wheel chair or medical bed.
- 3.Height range allows lifting from floor and high surfaces.
- 4.manual emergency stop for added safety
- 5.Adjustable knee pad unit for comfort at any height.
- 6.24V DC motor allows quiet, smooth operation
- 7.Warning beep for low battery charge and weight capacity levels
- 8.Lower base height maximizes stability and fits under most beds



Specifications

Weight Capacity	150kgs standard 200kgs optional
Overall Length	1100mm
Lifting range	390-1480/680-1280mm
Legs open width	840mm
Legs closed width	650mm
Front twin casters	2"
Overall height of legs Ground clearance	120mm
Rear break casters	4" with brakes
motor	220V,50HZ standard,8000N
Battery	24V Rechargeable sealed
Charger	Output: 24V DC
Warranty	Limited One (1) Year Warranty



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com

DG202

Battery-Powered Stand-up Lift



General Description :

WANROOEMED DG202 Battery-Powered Stand-up Lift is one ideal transport device to provide quick, safe assistance to those who have trouble standing. With multi-safety design, DG202 Battery-Powered Stand-up Lift ensures caregiver and resident security

Features :

- 1.Heavy-gauge steel construction
- 2.Leg open function to mostly close the wheel chair or medical bed.
- 3.Height range allows lifting from floor and high surfaces.
- 4.manual emergency stop for added safety
- 5.Adjustable knee pad unit for comfort at any height.
- 6.24V DC motor allows quiet, smooth operation
- 7.Warning beep for low battery charge and weight capacity levels
- 8.Lower base height maximizes stability and fits under most beds



Specifications

Weight Capacity	150kgs standard 200kgs optional
Overall Length	1200mm
Lifting range	760–1820mm
Legs open width	810mm
Legs closed width	640mm
Front twin casters	2"
Overall height of legs Ground clearance	120mm
Rear break casters	4" with brakes
Power	220V,50HZ standard
Battery	24V Rechargeable sealed
Charger	Output: 24V DC
Warranty	Limited One (1) Year Warranty



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanroomed.com www.wanroomed.com

Patient Lift Slings



DGA101

Economy All-purpose Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric
- 2.easy to install, adjust
- 3.excellent back support and protect
- 4.gray-blue
- 5.machine washable
- 6.loading 150KG

DGA406NE

Standard Grid Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric
- 2.easy to install, adjust
- 3.excellent back support and protect
- 4.gray-blue
- 5.machine washable
- 6.loading 150KG
- 7.It has large, medium and small size

DGA407NE

Standard Mesh Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric or green nylon mesh cloth
- 2.easy to install, adjust
- 3.excellent back and head support and protect
- 4.fit to buttocks area
- 5.leg pad
- 6.gray-blue
- 7.machine washable
- 8.loading 150KG
- 9.It has large, medium and small size

DGA400SF

Standing Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric
- 2.easy to connect and adjust
- 3.protect back support, it has fixed belt
- 4.suitable for standing type or multi-functional type patient lift
- 5.machine washable ; gray-blue ; 150KG ; loading 150KG
- 6.It has large, medium and small size



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanroomed.com www.wanroomed.com

Patient Lift Slings



DGA406N

Multi-functional Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric or green mesh fabric
- 2.easy to fit and adjust
- 3.not fit to buttocks and head area, excellent back protection
- 4.leg pad
- 5.it has fixed belt
- 6.gray-blue
- 7.loading 150KG
- 8.machine washable
- 9.It has large, medium and small size



DGA408N

Whole Body Protection Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric or green mesh fabric
- 2.easy to fit and adjust
- 3.whole body protection include head support and legs can be divided
- 4.suitable for the patient who needs whole body nursing
- 5.with toilet seat hole
- 6.leg pad
- 7.gray-blue
- 8.machine washable
- 9.loading: 150KG
- 10.It has large, medium and small size



DGA409N

Hole Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

Features :

- 1.green nylon knitted fabric or green nylon mesh cloth
- 2.easy to fit and adjust
- 3.excellent back and head support and protect
- 4.with toilet hole
- 5.gray-blue
- 6.loading 150KG
- 7.machine washable
- 8.It has large, medium and small size



**NOTICE: N type made of nylon knitting material, NE type made of nylon mesh material.
All available according to requirement.**



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com

DG201

Battery-Powered Patient Lift



General Description :

WANROOEMED Patient Lift is used for reducing the possibility of caregiver back injury and to ensure dignity in patient handling. These high-quality patient lifts are comfortable, reliable and integral tools for staff and caregivers. Using the patient lift for patient handling can help reduce staff injuries, save you money and reinforce safety as a priority in your health care facility

Features :

1. Leg base open easily by foot pedal
2. Height range allows lifting from floor and high surfaces.
3. manual emergency stop for added safety
4. Long, padded handles offer a plethora of grip choices.
5. Patient can turn 360° at any height during the lifting, to access a most comfort position.
6. The hoister handle design accords to dynamics, easy to push forward, backward, or turn.
7. Lower base height maximizes stability and fits under most beds
8. Warning beep for low battery charge and weight capacity levels



Specifications

Weight	43kgs
Weight Capacity	150kgs standard 200kgs optional
Overall Length	1100mm
Lifting range	740-1775mm
Legs open width	800mm
Legs closed width	585mm
Front twin casters	2"
Overall height of legs Ground clearance	120mm
Rear break casters	4" with brakes
Power	220V,50HZ standard
Battery	24V Rechargeable sealed
Charger	Output: 24V DC
Warranty	Limited One (1) Year Warranty
Package	1160*660*250mm, 144pcs/20GP, 324pcs/40GP, 360pcs/40HQ



Tel:0086-512-58810860 Fax:0086-512-58810663

Email: info@wanrooemed.com www.wanrooemed.com