


Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery



A-S106B

Electric Examination Table



江苏如翼医疗科技有限公司
JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.
Add: Room1311,1313, Guotai New World Plaza, No19,
East Renmin Road,Zhangjiagang, Jiangsu, China.
Tel:+86 512 58810663
Whatsapp/Mobile:+ 86-18962287398;15962399164
Email: info@wanrooemed.com
Website: www.rooemed.com www.wanrooemed.com



JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

A-S106B Electric Examination Table

A-S106B is suitable for all kinds of examinations including gynecological examination ,beautiful appearance,practical function and easy to operate.

Equipped with wired hand remote ,electric motor control tabletop elevation, gas spring control back plate.Manual control seat plate up and down.

Also equipped with recessed auxiliary plate, retractable foot pedal and waste basin.

The color of the mattress is optional. Foot switch and LED examination lamp are all optional.



DESCRIPTION ▶

- Electric motor control table top elevation.
- Integrated Waterproof mattress, adopt foaming technology, soft and comfortable.
- Back plate is gas spring control, strongly support ,easy to operate.
- Seat plate, auxiliary plate(waste basin in side), foot pedal all can manually control.
- Equipped with drawer, convenient for store.
- Equipped with hand remote, easy to operate table top elevation.
- Foot switch and examination light are optional.
- Mattress color are optional: blue, grey and black.



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

Model	Site	Length	Width	Tabletop Lowest Height	Tabletop Maximum height
A-S106B		1400mm	720mm	460mm	850mm

PARAMETER A-S106B Electric Examination Table

Movements		Model: A-S106B
Back Plate	up	70°
Seat Plate	up	7°
Auxiliary Plate Sliding Range		410mm
Foot Pedal Adjustment Angle Range		15° ~ 45°
Foot Pedal Outward Range		115°
Foot Pedal Sliding Range		350mm

OPTIONAL COLOR



STANDARD ACCESSORIES

Roll Stand	Waste Basin
1piece	1piece

Tabletop Lowest Height








Table Top Maximum Height




Foot Pedal Sliding Range




Foot Pedal Outward Range



Foot Pedal Adjustment Angle Range



Auxiliary Plate Sliding Range



A-S106B suitable for different kinds of examination, especially for gynecological.