

ALDR100B Hospital Birthing Bed

ALDR100B LDR bed (comfortable) nursing the lying-in woman **through all stages of parturition, including labor, delivery and recovery stage**. It can be easily changed into chair position to protect lying-in women from transfer from one bed to another. Electric controlled height adjustment, tilt adjustment and back section adjustment greatly enhance the work efficiency of the nurses.

Feature:

- Linak motor and electric control system.
- Centre brake and swivel system make the bed stable.
- Rotating side rails are easily operated.
- Minimum height enables patient get on the bed easily.
- Leg holder lockable, can turn forward & backward, left & right, height also adjustable.
- I.V. stand, height adjustable
- Head board, led board moveable.
- Cushion for newborn baby used during parturition.
- Stainless steel basin and flume moveable for washing.
- Grab handle for lying-in woman to exert her strength.
- Cushion ventilating, water proof, soft and three different colors available, i.e. rose, pink and blue.
- **Foot switch and hand control are optional.**

Technical Specifications:

BED SIZE	Overall Length:2050 mm
	Width: 980mm
	Height: 580-950mm(electric)
	Bed surface: 1815(L)*880mm(W)
FUNCTIONS	lowest≤580mm, height≥950mm (without the cushion)(electric)
	Trendelenburg: 0°-15°(electric)
	Reversed Trendelenburg: 0°-5°(electric)
	Backrest Adjustment: 0°-70°(electric)
	Leg section telescope: 550mm(manual)
Power supply	220V,50HZ, standard
Basic Accessories:	Leg support: 2 sets Clamp 1 pair Grab handle: 1 pair Foot rest: 1 pair Cushion: 1 set Filth Basin: 1 set IV pole: 1 set hand controller: 1pc
Package	Wooden package(1pc/case): 187*112*110cm, Net weight: 200kg,Gross weight: 290kg

Accessories Trolley(optional)



ALDR100BM Hydraulic Hospital Birthing Bed

ALDR100BM Hydraulic Hospital Birthing Bed is economic birthing bed for the lying-in woman **through all stages of parturition, including labor, delivery and recovery stage**. It can be easily changed into chair position to protect lying-in women from transfer from one bed to another. Hydraulic controlled height adjustment, manual tilt adjustment and back section adjustment greatly enhance the work efficiency of the nurses.

Feature:

- USA oil pump supply the power source for height adjustment
- tilt adjustment and back section
- Centre brake and swivel system make the bed stable.
- Rotating side rails are easily operated.
- Minimum height enables patient get on the bed easily.
- Leg holder lockable, can turn forward & backward, left & right, height also adjustable.
- I.V. stand, height adjustable
- Head board, foot board moveable.
- Cushion for newborn baby used during parturition.
- Stainless steel basin and flume moveable for washing.
- Grab handle for lying-in woman to exert her strength.
- Cushion ventilating, water proof, soft and three different colors available, i.e. rose, pink and blue.

Technical Specifications:

- Bed dimension including bed end and bumper wheel: Length 2050mm, Width 980mm
- Bed dimension: Length 1815mm, Width 880mm

- Bed height: lowest \leq 580mm, height \geq 950mm (without the cushion)(hydraulic)
- Back section turning: \geq 70° (gas spring)
- Reversed Trendelenburg: \geq 5°(gas spring)
- Trendelenburg \geq 12°(gas spring)

Leg section telescope: 550mm(manual) Castor dia.: 125mm

Standard accessories:

Leg holder 1 pair

Clamp 1 pair

Grab handle 1 pair

Foot rest 1 pair

Basin 1 pc

Water proof cushion 1 set

I.V. stand 1 pc

Cushion for newborn baby: 1 set hand controller: 1pc

Cushion: 1 set (includes back cushion,Seat cushion and assist platform cushion)

Accessories Trolley(optional)

