



# TECHARTMED®

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

## ET200 Ophthalmological operating table + electric doctor chair

The table is designed and made specially according to ophthalmologists' advices and can apply to various operations for ophthalmology.  
Tabletop lift adopts electric motorized type ,easy to operate.  
All cover made of 304 stainless steel,easy to clean.  
Nice appearance and stable performance can be regarded as an ideal ophthalmological operating table.  
Equipped with electric chair,makes the doctor more comfortable and convenient in operation.



Your Ideal Operating  
Table Supplier


Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery



# TECHARTMED®

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery



## ET200

Ophthalmological operating table  
+ electric doctor chair



# TECHARTMED®

NINGBO TECHART MEDICAL EQUIPMENT CO.,LTD

Add: No.99 South Hongtang Road Ningbo China 315033  
Tel: +86-574-83076915  
Fax: +86-574-87137656  
E-mail:sales01@techartmed.com sales02@techartmed.com sales03@techartmed.com  
Web: www.techartmed.com

NINGBO TECHART MEDICAL EQUIPMENT CO., LTD





# TECHARTMED®

## Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

## Production details

### Ophthalmological operating table

- Foot switch: several groups of different functions is convenient for operation by doctor.
- Back plate and leg plate: adopted the foot switch that enables flexible and smooth adjustments.
- Head plate: is easy to control by mechanical gear.



### Electric chair

- Adjustable armrest support with rotation 180°
- Adjustable backrest in height.
- Chair elevation(up/down): two buttons on the chair base,the doctor can be adjusted by foot.
- Double brakes system on the wheels.



## Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery



# TECHARTMED®

## Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

Model	Site	Length	Width	Tabletop lowest height	Elevation range
ET200		≥2080mm	≥550mm	≤500mm	≥200mm

PARAMETER ET200 ophthalmological operating table + electric doctor chair

Movements		Model:ET200
Head plate	Up	≥45°
	Down	≥90°
Leg plate	Up	≥30°
	Down	≥30°
Back plate	Up	≥45°
	Down	≥20°



## OPTIONAL COLOR



### ELECTRIC CHAIR

Model	Site	Chair top diameter	Minimum height	Elevation range
Electric chair		420mm	460mm (chair top)	140mm

