

*support you care you*



**WANROOEMED**

WELL SOLD FOR 40 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL SURGERY

**AOT400**

**AOT400A**



**AOT400C**

**support you care you**

## AOT400

WELL SOLD FOR 40 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL SURGERY



### Introduction

Endearing appearance, various colors of mattress available., optional pedal switch electrically and easily operate.

It has overcome the slowness matter of similar products.

### Main Technical Data

Parameters may have reasonable error.

AOT400A	AOT400C
Length and width for tabletop : 1820/520mm	Length and width for tabletop : 1780/600mm
Elevation(down/up) : 680/1090mm	Elevation(down/up) : 680/880mm
Trendelenburg/anti-trendelenburg : 14°/17°	Trendelenburg/anti-trendelenburg : 8°/20°
Back plate(up/down):60°/10°	Back plate(up/down):45°/7°
Power : 220V , 50Hz,1.0kW	Power : 220V , 50Hz,1.0Kw

*support you care you*



**WANROOEMED**

WELL SOLD FOR 40 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL SURGERY

**AOT400B、AOT400BM**

**AOT400B**



**AOT400BM**



*support you care you*

## AOT400B, AOT400BM

WELL SOLD FOR 40 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL  
SURGERY



### Introduction

400B adopted both remote control and pedal switch, diversified operation.

The cover for base and column is made of 304 stainless steel, good performance on anti-rust and easy to clean. 400BM with endearing appearance, seamless mattress, mechanically hydraulic control. Backplate is adjusted by germany made gas spring.

### Main Technical Data

Parameters may have reasonable error.

AOT400B	AOT400BM
Length and width for tabletop : 1780/600mm	Length and width for tabletop : 1780/600mm
Elevation(down/up) : 680/880mm	Elevation(down/up) : 680/880mm
Trendelenburg/anti-trendelenburg : 8°/20°	
Back plate(up/down):45°/7°	Back plate(up/down):45°/10°
Power : 220V,50Hz,1.0kW	