

## Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

### **AOT500** Electric operating table

Developments made in the conjunction with doctors: this is one of the basic principles which has made WANROOEMED the world market leader for operating tables. Close cooperation and years of experience have resulted in tried and tested products which facilitate the work in the operating room and help to increase the safety and comfort of the patient.

Best example: this table can withstand even extreme loads in any position. The structure enables it to be adapted to match the size of the patient and in no time at all it can be equipped for any surgical discipline. A longitudinal shift of up to 350mm makes the operating table even more versatile.



## Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

# WANROOEMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

Electric operating table

## AOT500



WANROOE MACHINERY CO., LTD.

# WANROOEMED

Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery


## AOT500 Electric operating table



### Product Description

Engineered for increasingly complex surgical scenes, AOT500 comes with intelligent automatic control system satisfying various demands of general surgeries. Extra wide tabletop, extreme loading system, long horizontal sliding that can be suitable for both X-ray and C-arm. AOT500 is an ideal electric operating table.

- Longitudinal slide up to 350mm for free access to C arm
- Extreme weight load capacity of 250kgs with highest safety and stability
- Memory pad with antistatic, waterproof design
- Standard override panel control makes the control safer
- Battery power for 50-80 operations



### PARAMETER AOT500 Electric operating table

Movements		Model: AOT500
Tabletop trendelenburg	Reverse Trendelenburg	≥ 25°
	Trendelenburg	≥ 25°
Lateral Tilt (Left And Right)		≥ 15°
Head plate	Up	≥ 40°
	Down	≥ 90°
	Flexible	≥ 40°
Back plate	Up	≥ 80°
	Down	≥ 20°
	Up	≥ 15°
Leg plate	Down	≥ 90°
	Outward	≥ 90°
	Kidney Bridge Up	110mm
Horizontal sliding		≥ 350mm

### Details

1. Longitudinal shift: the longitudinal shift of up to 350MM enables optimum access with the C-arm without needing to reposition the patient or to reverse the positioning of the patient. Each selected fluoroscopy point may easily be moved in front of or behind the column and ensures comfortable access to the patient from initiation to surgery.
2. It can easily be moved to any desired position as four castors.
3. Leg plate: imported gas spring for easy controls.
4. Two joint head plate.
5. Built-in kidney bridge.
6. Double control system.



# Operating Table

Well Sold For 40 Years Worldwide  
Meet The Basic Demands Of General Surgery

Your Ideal Operating Table Supplier

Model	Site	Length	Width	Tabletop lowest height	Elevation range
AOT500		≥ 2100mm	≥ 550mm	≤ 710mm	300mm

### STANDARD ACCESSORIES

Anaesthesia screen	1 piece	Body support	1 pair	Leg support	1 pair	Hand remote	1 piece
Kidney bridge handle	1 piece	Arm support	1 pair	Foot support	1 pair	Mattress	1 set
Shoulder support	1 pair	Long fixing clamp	1 pair	Fixing clamp	8 pieces	Power Line	1 piece

