

A good operating table for your OR room

## AOT700 Electric-hydraulic Operating Table



Wanrooe Machinery Co., Ltd

Add: Jinfeng Town, Zhangjiagang, Jiangsu, China

Tel: 0086-512-58810860

Fax: 0086-512-5881063

Website: [www.wanrooemed.com](http://www.wanrooemed.com)

Email: [info@wanrooemed.com](mailto:info@wanrooemed.com),

[info@pallmanngroup.com](mailto:info@pallmanngroup.com)

**WANROOEMED, MAKE THE QUALITY OPERATING TABLE**

## AOT700 Electro-hydraulic Operating Table

### Feature:

Silent Electric-hydraulic system, 304 stainless steel base and cover, all accessories made of quality 304 stainless steel

C-Arm compatible table top

225 KG (500 LB) lift capacity.

Horizontal slide motion to create a greater field of view for C-arm use.

Kidney Bridge elevation

Mobile anti-static castors, controlled by hand controller

"Flex" positioning

Battery backup

Optional: Double remote control with auto-reset / level function







# AOT700 Electric-Hydraulic Operation Table





# AOT700 Electric-Hydraulic Operation Table



选配1005牵引架胫骨骨钉  
接入手术  
Optional 1005 Tibia  
bone-nailing position



选配1005牵引架平躺位  
Optional traction 1005 supine position



选配1005牵引架侧卧位  
Optional 1005 lateral lying position



选配1006牵引架  
Optional traction 1006



选配1005牵引架胫骨手术体位  
Optional 1005 Tibia extension position



1005牵引架接口  
Joint for 1005 lateral lying position



选配脊椎手术架  
Optional spondyle operating  
holder



## AOT700 Electric-Hydraulic Operation Table



选配侧卧搁手板  
Operation lateral arm support



选配妇科污物盆  
Operation basin for gynecology



胸、肾、胆囊侧卧手术台  
Lateral lying position for chest,  
kidney, cholecyst operation



肠镜手术位  
Anorectum operation position



耳鼻喉手术位  
E. N. T operation position



胸腹手术位  
Chest and Abdomen operation  
position

# AOT700 Electric-Hydraulic Operation Table



X-ray双层面板 Operation X-ray Board



台面平移 Tabletop headward sliding/toward leg sliding



台面升降 Tabletop elevation up/down



8CM记忆床垫  
8cm thickness  
memory foam mattress



内置腰桥  
Built-in kidney bridge



选配双控遥控(含一键复位)  
Optional second remote control  
with reposition function



头板调整(上/下)  
Head plate adjustment  
(up/down)



腿板调整(上/下)  
Leg plate adjustment  
(up/down)



腿板调整(向外)  
Leg plate adjustment  
(outward)



左/右倾 Lateral tilt left/right



前/后倾 Trendelenburg/anti-trendelenburg



台面弯曲 Tabletop flexural position



刹车 Brake system lock/release

## 主要规格参数 Main specification

注：表列参数与实物可能会有合理误差

项目	Item	内容 Contents
台面长度及宽度	Length and width for tabletop	1960mm×500mm
台面最低及最高高度	Elevation(up/down)	690mm/1100mm
台面前倾后倾最大角度	Trendelenburg/Reverse Trendelenburg	30°/30°
台面左右倾最大角度	Lateral Tilt(left/right)	18°/18°
头板转折	Head rest(up/down)	40°/90°
腿板转折最大角度	Foot Plate(up/down&horizontal positioning)	15°/90° 90°
背板转折最大角度	Back Plate(up/down)	90°/40°
台面平移	Longitudinal Sliding	350mm
腰桥升距	Wrist Plate up	110mm
电源电压、电源频率、电源容量	Power	~220V, 50Hz, 1.0kW

## 标准配件 Standard accessories

麻醉屏架	1根	肩托	1对	搁手板	1对	支身架	1对
托腿架	1对	长滑块	1对	足板	1对	手控器	1个
床垫	1套	内置腰桥摇手柄	1个	电源线	1根	液压油	1壶

