


## F-PS90L Plasma Sterilizer



### Features:

1. 1.7 inch LCD(liquid crystal display) and Android operating system;
2. Finger print identification and no upper limit to the number of operators;
3. Support wired LAN,WLAN and 3G;
4. 2 USB connectors,which can export data record;
5. Text message reminder;
6. Video broadcast;
7. Cloud service,backup and recovery record data;
8. Real time curve:the curve shows the sterilization status for easy monitoring;
9. Online infection control knowledge,user manual and maintenance instructions are available;
10. update online medical industry news;
11. Update application online and offline;
12. Export record data as excel file;
13. Leybold vacuum pump from Germany.



*to order*

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com)      TEL: 0086-512-58810860      FAX: 086-512-58810663

Type	F-PS90L
Overall Dimension	H1725*W820*D890(mm)
Door Type	Single automatic vertical sliding door
Capacity	120L(cylinder chamber)
Chamber Size	Φ460*770(mm)
Sterilization Mode	Quick/standard/enhanced cycle
Sterilization Temp	50°C--60°C
Power Rating	220V/50Hz,2.7Kw
Control system	Android system,intelligent management
Cassette	12 capsules for each cassette which can be used for 6 cycles

WANROOEMED

to order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663