

豪华型午休椅

Luxurious Chairbed

专利号: ZL 2012103336053 (发明专利)

ZL 2012304332830

ZL 2012204600665

ZL 201220460028X

好午休健康添活力
午休好工作更给力



您想要的就是我们正在做的
What you want is right what we are doing now

型号:BHC001M
名称:豪华午休椅(陪护椅)
Luxurious Noon break Chair
规格(cm):83×89×45(靠背高:100)
展开(cm):198×89×45
净重:70kg

说明:面料为环保型PVC压延仿皮, 钢管喷塑框架, 内脚为高密度板包皮脚板, 隐藏式静音脚轮, 高档实木扶手, 可调节靠枕, 椅子内部有伸缩弹簧, 收放自如。

Description:The Fabric is pressed PVC artificial leather, steel frame coated, foot board made of high density board and leather, hidden mute casters, wooden armrest, the pillow can adjust, extension spring inside, easy to close and open.



Tel:0086-512-58810860

www.wanrooemed.com

Email: info@wanrooemed.com

豪华型午休椅

Luxurious Chairbed

专利号: ZL 2012103336053 (发明专利)

ZL 2012304332830

ZL 2012204600665

ZL 201220460028X

好午休健康添活力
午休好工作更给力



型号:BHC001D
名称:豪华午休椅 (陪护椅)
Luxurious Noon break Chair
规格(cm):80×75×46(靠背高:100)
展开(cm):198×75×50
净重:48kg

说明:超纤维环保人造革面料, 钢管喷塑框架, 高密度板包皮脚板, 隐藏式静音脚轮, 钢木结构扶手, 可调节靠枕, 椅子内部有伸缩弹簧, 收放自如。

Description:The Fabric is microfiber leather, steel frame coated, foot board made of high density board and leather, hidden mute casters, wooden and steel armrest, the pillow can adjust, extension spring inside, easy to close and open.



隐藏式脚轮



靠枕0-18° 调节

您想要的就是我们正在做的
What you want is right what we are doing now



Tel:0086-512-58810860

www.wanrooemed.com

Email: info@wanrooemed.com

豪华型午休椅

Luxurious Chairbed

专利号: ZL 2012204821704

创造更多
办公空间



型号:BHC001C
名称:豪华午休椅(陪护椅)
Luxurious Accompanier's Chair
规格(cm):78×65×45(靠背高:106)
展开(cm):206×65×45
包装尺寸(cm):130×66×38
净重:38kg



说明:铝合金扶手, 静音脚轮, 面料为环保型PVC压延仿皮, 钢管喷塑, 适合高级病房及办公室。
Description:aluminum alloy armrest, mute wheels,the fabric is pressed PVC artificial leather,for luxurious sickroom and office.



您想要的就是我们正在做的
What you want is right what we are doing now



Tel:0086-512-58810860 www.wanrooemed.com Email: info@wanrooemed.com



型号: BHC001E
 名称: 豪华陪护椅
 Luxurious Accompanier's Chair
 规格(cm): 88×70×62(靠背高:100)
 展开(cm): 200×70×46
 净重: 26kg
 说明: 面料为环保型PVC压延仿皮, 静音脚轮, 钢管喷塑框架, 内有气弹簧, 收放简单, 椅面造型设计。
 Description: the fabric is pressed PVC artificial leather, mute wheels, steel frame coated, gas spring inside, easy to close and open, chair surface modelling design.



型号: BHC001K
 名称: 豪华陪护椅
 Luxurious Accompanier's Chair
 规格(cm): 80×71×43(靠背高: 98)
 展开(cm): 204×71×43
 净重: 36kg
 说明: 面料为PVC压延仿皮, 实木扶手, 钢管喷塑。
 Description: the fabric is pressed PVC artificial leather, wooden armrest, steel frame coated.

型号: BHC001J
 名称: 豪华陪护椅
 Luxurious Accompanier's Chair
 规格(cm): 78×75×45(靠背高:106)
 展开(cm): 206×75×45
 包装尺寸(cm): 130×66×38
 净重: 36kg
 说明: 面料为PVC压延仿皮, 聚氨酯结皮发泡扶手, 钢管喷塑。
 Description: the fabric is pressed PVC artificial leather, Polyurethane foam armrest, steel frame coated.



型号: BHC001L
 名称: 豪华陪护椅(加宽型)
 Luxurious Accompanier's Chair
 规格(cm): 80×102×43 (靠背高: 98)
 展开(cm): 208×102×43
 净重: 55kg
 说明: 面料为亚麻布, 钢管喷塑。
 Description: the fabric is linen, steel frame coated.



您想要的就是我们正在做的
 What you want is right what we are doing now



型号:BHC001F
 名称:新型陪护椅
 New Accompanier's Chair
 规格(cm):66×68×41(靠背高:93)
 展开(cm):186×68×40
 净重:32kg
 说明:面料为高档PVC压延仿皮, 钢管喷塑, 内有伸缩弹簧, 收放自如。
 Description:the fabric is pressed PVC artificial leather,extension spring inside,easy to close and open.



型号:BHC001S
 名称:不锈钢陪护椅
 Stainless Steel Accompanier's Chair
 规格(cm):65×62×42
 展开(cm):190×62×58
 包装尺寸(cm):118×62×22
 净重:18kg
 说明:面料为PVC压延仿皮。
 Description:the fabric is pressed PVC artificial leather.



型号:BHC001
 名称:陪护椅
 Accompanier's Chair
 规格(cm):65×62×42
 展开(cm):190×62×58
 包装尺寸(cm):118×62×22
 净重:22.5kg
 说明:面料为PVC压延仿皮, 钢管喷塑。
 Description:the fabric is pressed PVC artificial leather.

可选配件: 锁

型号:BHC001B
 名称:双折陪护床
 Double Fold Accompanier's Bed
 规格(cm):44×55×98
 展开(cm):182×55×42
 包装尺寸(cm):58×46×100
 净重:20kg
 说明:面料为沙发用棉布, 钢管喷塑。
 Description:the fabric are the cotton fabrics sofas.



型号:BHC001A
 名称:可藏式陪护椅
 Hiding Accompanier's Chair
 规格(cm):74×62×35
 展开(cm):190×62×20
 包装尺寸(cm):74×66×38
 净重:21kg
 说明:面料为PVC压延仿皮, 钢管喷塑。
 Description:the fabric is pressed PVC artificial leather.

您想要的就是我们正在做的
 What you want is right what we are doing now



型号:BHC001H
名称:实木陪护椅
Wooden Accompanier's Chair
规格(cm):65×62×42
展开(cm):190×62×58
包装尺寸(cm):118×62×22
净重:46kg
说明:面料为橡木皮, 钢管喷塑。
Description:the fabric is oak board.

型号:BHC003E
名称:豪华气动输液椅
Luxurious Gas spring Transfusion-Chair
规格(cm):80×80×46 (靠背高: 113)
净重:28kg
说明:面料为环保型PVC压延仿皮, 钢管喷塑, 内有气弹簧, 任意角度调节。
Description:the fabric is pressed PVC artificial leather,gas spring inside any angles adjustable.



型号:BHC001G
名称:陪护床(加床)
New Accompanier's Bed
规格(cm):75×24×98
展开(cm):186×75×30
包装尺寸(cm):70×24×92
净重:18.5kg
说明:绒布面料, 钢管喷塑。
Description:the fabric is lint or leather,steel frame coated.



型号:BHC003A BHC003AE
名称:豪华型输液椅
Luxurious Transfusion-Chair
规格(cm):90×75×110
包装尺寸(cm):100×60×44
净重:31kg
说明:面料为环保型PVC压延仿皮, 钢管喷塑, 有三种角度可以调节, 带旋转餐桌。
Description:the fabric is pressed PVC artificial leather,3 angles adjustable,with rotatable table.

型号: S-7E
名称: 豪华电动输液椅
Luxurious Electric Transfusion-Chair

您想要的就是我们正在做的
What you want is right what we are doing now



型号:BHC003
名称:常规输液椅
Transfusion-Chair
规格(cm):77×68×108
包装尺寸(cm):92×62×35
净重:21.5kg
说明:面料为PVC压延仿皮, 钢管喷塑, 有三种角度可以调节。
Description:the fabric is pressed PVC artificial leather,3 angles adjustable.



型号:BHC003B
名称:固定式输液椅
Transfusion-chair
规格(cm):77×68×108
包装尺寸(cm):82×62×30
净重:19kg
说明:面料为PVC压延仿皮, 钢管喷塑, 不可调节角度。
Description:the fabric is pressed PVC artificial leather,not adjustable.



型号:BHC003C
名称:不锈钢输液椅(抽拉式)
Stanless Steel Transfusion-Chair
规格(cm):77×68×108
包装尺寸(cm):92×62×35
净重:20kg
说明:面料为PVC压延仿皮, 二种角度可调节。
Description:the fabric is pressed PVC artificial leather, 2 angles adjustable



输液室效果图



型号:BHC003S
名称:不锈钢输液椅
Stanless Steel Transfusion-Chair
规格(cm):77×68×108
包装尺寸(cm):92×62×35
净重:22kg
说明:面料为PVC压延仿皮, 带餐桌(选配), 三种角度可调节。
Description:the fabric is pressed PVC artificial leather,with table, 3 angles adjustable.



型号:BHC003D
名称:新型输液椅
New Transfusion-Chair
规格(cm):82×70×102
包装尺寸(cm):90×56×30
净重:25kg
说明:面料为环保型PVC压延仿皮, 钢管喷塑, 角度不可调节。
Description:the fabric is pressed PVC artificial leather,not adjustable.



型号:BHC003DC
名称:儿童输液椅

您想要的就是我们正在做的
What you want is right what we are doing now