


DB863-I Electric Portable Thermal Jade Massage Bed




1. 8+1 infrared heated jade rollers for lumbar.
2. Synchronous belt to avoid track slip.
3. Lower section 4 kneading rollers.
4. 5 balls jade projector to massage all sections of your body.
5. Six auto modes plus one manual mode, 40 minutes auto timer.
6. Temperature is adjustable from 40-70 centigrade degree.
7. Liquid Crystal Display remote controller.

Leather color: optional to your requirement.

Carton size: 113x57x22cm

Packaging: 1pc/ctn CBM: 0.142 N.W./G.W.: 31/33kg

20' QTY.: 200pcs 40' QTY.: 400pcs 40' HC QTY.: 470pcs



Details:

Jade Massage Bed



5 Balls Hand massage head



LCD hand controller

Optional: High quality steel frame(large weight,easy to move)



Jade Massage Bed

DB863TJT-YZP Portable Jade Massage Bed



- 1.Upper section 5 wide jade rollers(for neck/head and lumbar massage).S Shaped Patented Drive Rail in upper back will follow the contour of your spine!
- 2.Whole body carbon fiber heating(FIR)for ultimate therapy.
3. 5 + 11 jade projectors to massage all parts of your body.
4. Six auto modes plus one manual mode,40 minutes auto timer.
- 5.Temperature is adjustable from 40-70 centigrade degree.
- 6.Liquid Crystal Display remote controller.
- 7.Folded design for more room space.
- 8.Metal stand with universal wheel.
9. Large Bed Size will work for anyone up to 6.5' feet tall! Solid Heavy Duty metal construction allows any user up to 350 pounds.

NOTE: MP3 optional

Package size: 138x76x48cm packing: 1pc/1ctn N.W./G.W. 63/67kg CBM 0.5 20'QTY. 56PCS 40'QTY.
112PCS 40'HC QTY. 136PCS

Folded:



Jade Massage Bed



Model DB863JT-R Massage Bed



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupressure, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

S-shaped driving orbit, aluminum rail can be selected

LCD Display provides six auto-working modes

Equipped with 11pcs jade balls+5pcs jade balls

Five far-infrared jade rollers in the upper body

Legs kneading function

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ

Voltage:120W

Size: 198* 65* 49cm

Package: 105*67.5*22.5(bedstead)

125*66.5*31.5(bed)


Gross weight: 75kgs



Jade Massage Bed



Model DB100 Massage Bed



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupressure, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

S-shaped driving orbit, aluminum rail can be selected

Synchronous belt transmission, more durable

Whole body with nine jade rollers, upper 5 rollers, lower 4 rollers

LCD controller shows six auto-working modes

Equipped with 11 pcs jade balls + 5 pcs jade balls

Carbon fiber heating for the upper body

Optional: LCD TV which can play MP3, MP4 ect.

Technical Specifications:

Power supply: AC220V, 50HZ or AC110V, 60HZ

Voltage: 120W Size: 198* 65* 49cm

Package Details:

2 Cartons/pc, 0.42CBM/pc, 56pcs for 20 ft, 114pcs for 40 GP, 150pcs for 40HQ

Bed: 125X66.5X31.5 cm G.W.:50Kg N.W.:46Kg

Bedstead :105X67.5X22.5cm G.W.:42Kg N.W.:38Kg

Jade Massage Bed

5 Balls Hand massage head

5 balls hand massage head with the five natural jade heating function composition. According to a certain extent word order, crescent-type can be a good snapping the body. Jade massage heads can be fever, the temperature can be adjusted according to your needs.

11 balls hand massage head

11 balls hand massage head with heating function by the 11 natural jade components. According to a certain extent was arranged in an array, can be a good snapping the body. Jade massage heads can be fever, the temperature can be adjusted according to your needs.



back massage head

back massage head composed of 4 back massage rollers, 1 neck massage roller and 1 jade cap. The massage heads is heated by 11 halogen lamps to a temperature between 30 °C-70 °C adjustable. Massage head massage program under the control of the human spinal manipulative correction massage, warm round jade massage can promote blood circulation, relieve pain, fatigue disperse your day



Massage Rollers for the leg



4 massage heads and one jade cap, heated by 9 Halogen lights, temperature adjustable from 30°C to 70°C

LCD hand controller

Hand controller using a large-screen LCD, graphic display is clear and intuitive



Jade Massage Bed



Model DB101 Massage Bed(all-in-one machine)



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupuncture, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

Designed with straight and **aluminum rail**

Driving by synchronous belt, more durable

Whole body massage with nine jade rollers, upper 5 rollers, lower 4 rollers and carbon fiber heating

LCD controller shows six auto-working modes

Equipped with 11pcs jade balls + 5pcs jade balls

MP3 player

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ

Voltage:120W v Size: 198* 65* 49cm

Package Details:

2 Cartons/pc, 0.46CBM/pc,

Bed:198X66.5X24 cm G.W.:49.5Kg N.W:45.5Kg

Bedstead:105X67.5X22.5cm G.W.34Kg N.W.30Kg

50pcs for 20 ft, 105pcs for 40 GP, 135pcs for 40HQ

Jade Massage Bed

5 Balls Hand massage head

5 balls hand massage head with the five natural jade heating function composition. According to a certain extent word order, crescent-type can be a good snapping the body. Jade massage heads can be fever the temperature can be adjusted according to your needs.



11 balls hand massage head

11 balls hand massage head with heating function by the 11 natural jade components. According to a certain extent was arranged in an array, can be a good snapping the body. Jade massage heads can be fever, the temperature can be adjusted according to your needs.

back massage head



Aluminum alloy straight rails, more durable.

back massage head composed of 4 back massage rollers, 1 neck massage roller and 1 jade cap. The massage heads is heated by 11 halogen lamps to a temperature between 30 °C-70 °C adjustable. Massage head massage program under the control of the human spinal manipulative correction massage, warm round jade massage can promote blood circulation, relieve pain, fatigue disperse your day


Massage Rollers for the leg 4 massage heads and one jade cap, heated by 9 Halogen lights, temperature adjustable from 30°C to 70°C

LCD hand controller

Hand controller using a large-screen LCD, graphic display is clear and intuitive



MP3 music



Model DB102 Massage Bed



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupressure, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

S-shaped driving orbit, more comfortable

LCD controller, equipped with 11 pcs jade balls + 5 pcs jade balls

Five far-infrared jade rollers for the back with neck massage

Whole body heating through the carbon fiber and tourmaline

Backrest could be ascended and descended freely

Technical Specifications:

Power supply: AC220V, 50HZ or AC110V, 60HZ

Voltage: 120W Size: 198* 65* 49cm

Package Details:

2 Cartons/pc, 0.46CBM/pc,

Bed: 125X66.5X31.5 cm G.W.: 42.5Kg N.W.: 38.5Kg

Bestead: 129X66.5X23.5cm G.W.: 54Kg N.W.: 50Kg

Jade Massage Bed

51pcs for 20 ft, 105pcs for 40 GP, 133pcs for 40HQ

5 Balls Hand massage head

5 balls hand massage head with the five natural jade heating function composition. According to a certain extent word order, crescent-type can be a good snapping the body. Jade massage heads can be fever the temperature can be adjusted according to your needs.



11 balls hand massage head

11 balls hand massage head with heating function by the 11 natural jade components. According to a certain extent was arranged in an array, can be a good snapping the body. Jade massage heads can be fever, the temperature can be adjusted according to your needs.

back massage head

back massage head composed of 4 back massage rollers, 1 neck massage roller and 1 jade cap. The massage heads is heated by 11 halogen lamps to a temperature between 30 °C-70 °C adjustable. Massage head massage program under the control of the human spinal manipulative correction massage, warm round jade massage can promote blood circulation, relieve pain, fatigue disperse your day




Legs carbon fiber heating

Leg using carbon fiber heating with heat, long life, far infrared emission of the closest to human. Carbon fiber heating plate temperature between 30 °C to 70 °C adjustable.



Germanite

Germanium stone is umbrella name of the nature of the element germanium which contain a class of ore. Germanium is a light-colored metal, semiconductor materials, which atom has 32 electrons, while the expansion of the role of electronic engineering in the body, but also increase oxygen to the cells through the purification. Germanium restrains radiation-induced damage, thereby reducing damage and restoring the injured cells. Blood cells to increase oxygen making the blood clean



LCD hand controller

Hand controller using a large-screen LCD, graphic display is clear and intuitive



Model DB103 Massage Bed



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupuncture, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

S-shaped driving orbit, aluminum rail can be selected

Digital Display provides five auto-working modes

Equipped with 11pcs jade balls+5pcs jade balls

Five far-infrared jade rollers in the upper body

Carbon fiber heating for the leg

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ

Voltage:120W Size: 198* 65* 49cm

Package Details:

2 Cartons/pc, 0.35CBM/pc,

G.W.:60kgs, N.W.:56kgs

65pcs for 20 ft, 139pcs for 40 GP, 175pcs for 40HQ



Model DB104 Massage Bed



This product is our company based on the principles of Chinese massage and meridian theory, the use of natural jade material, developed from the application of modern high-tech. She is the integration of chiropractic, acupressure, massage, far infrared light waves, warm high-tech health products. It is a good choice for spine health care, especially for the elderly care massage products.

The massage table comes with a unique design hand massage ball, can be easily used in various parts of the body. The all-natural Xiuyan ball massage stone, heating inside the body which can be adjusted according to temperature, giving you the warmest care!

S-shaped working orbit, aluminum rail can be selected

Nine jade rollers inside for the whole body

LCD controller shows six auto-working modes

Equipped with 11pcs jade balls+5pcs jade balls

Carbon fiber heating for the upper body

Backrest with an Auto-Lifter ascend and descend freely

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ Voltage:120W Size: 198* 65* 49cm

Package Details:

2 Cartons/pc, 0.46CBM/pc, G.W.:98kgs, N.W.:94kgs

125*66.5*31.5cm(bed) 129*66.5*23.5cm(bedstead)

51pcs for 20 ft, 105pcs for 40 GP, 135pcs for 40HQ

Jade Massage Bed



Model DB104 Massage Bed



www.wanrooemed.com

- S-shaped driving orbit, aluminum rail can be selected
- Synchronous belt transmission, more durable
- Whole body with nine jade rollers, upper 5 rollers, lower 4 rollers
- LCD controller shows six auto-working modes
- Equipped with 11pcs jade balls+5pcs jade balls
- Carbon fiber heating for the upper body
- Backrest could be ascended and descended freely

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ

Voltage:210W

Package Details:

Bed 125*66.5*31.5cm, G.W, 50kgs N.W.: 46Kgs

Bed stand: 129*66.5*23.5cm, G.W, 54kgs N.W.: 50 kgs

51PCS/20GP, 105PCS/40GP, 133PCS/40HQ



www.wanrooemed.com

TEL: 0086-512-58810860

FAX: 0086-512-58810663

E-mail:info@wanrooemed.com

[Http://www.wanrooemed.com](http://www.wanrooemed.com)

Model DB104 Music Jade Massage Bed

智能温玉保健床 Jade Massage Bed



S-shaped working orbit, aluminum rail can be selected
Nine jade rollers inside for the whole body
LCD controller shows six auto-working modes
Equipped with 11pcs jade balls+5pcs jade balls
Carbon fiber heating for the upper body
Backrest with an Auto-Lifter ascend and descend freely

With MP3 Function

Technical Specifications:

Power supply: AC220V,50HZ or AC110V,60HZ Voltage:210W

Package Details:

Bed 125*66.5*31.5cm, G.W, 50kgs N.W.: 46Kgs

Bed stand: 129*66.5*23.5cm, G.W, 54kgs N.W.: 50 kgs

51PCS/20GP, 105PCS/40GP, 133PCS/40HQ

DB300 Korea Infrared Spine Protatable V3 Jade Roller Massage bed



Functions:

- 1.S-shaped walking orbit ,aluminum rail can be selected
2. Control panel is show by LCD,shows six auto-working modes
- 3.Equipped with 11 pcs jade balls + 5 pcs jade balls
- 4.four BALL jade rollers in the upper body,the massage roller can adjust the height can be adjusted .
- 5.The default system of temperature is 40°C.It can be ascended to 70°C according to your requirements.
- 6.Different colors of covers are available due to customers.
- 7.Body of bed weight:35kg;Leg of bed weight:40kg
Body of bed size:125*67*19cm;Leg of bed size:143*72*26cm

Package Size:	127X68X19cm 143x70x26cm
Qty. per Ctn.:	1 piece / carton
CBM for one bed:	0.424
20', 40', 40'H Qty.:	66, 132, 162 pieces
Rated Voltage:	100~240V, 50~60HZ
Rated Power:	150w
Certifications:	Medical CE, RoHS

DB400 message chair

Jade Massage Bed



- 1, liquid crystal display, easy to understand, easy to operate,
- 2, zero gravity setting,
- 3, L-track set more in line with human body structure, massage to the lower body,
- 4, increase jade massage, massage area larger,
5. Temperature 30 to 70 degrees adjustable, to meet the needs of different consumers,
- 6, the combination of leg airbags and wheels, to achieve a full range of massage function.

DB500 message mat



锺石按摩垫

排毒养颜，排除体内残留的重金属，使细胞活性化。

消除疲劳提高睡眠质量。

疏通经络，缓解各种人体神经性疼痛、缓解腰背及关节疼痛。

扩张微血管，促进血液循环，提高人体免疫力。

Detox beauty, remove the body of residual heavy metals, so that cell activation.

Eliminate fatigue and improve sleep quality

Clear the meridians, relieve a variety of human neuropathic pain, relieve back and joint pain.

Dilate capillaries, promote blood circulation, improve human immunity.

DB600 Vibration Massage Bed



振动床垫

Vibration Massage Bed

- 1、微电脑控制，自动定时。
 - 2、八种自动模式，三种强度随意选择。
 - 3、过流、过热保护。
 - 4、低压设计，安全舒适。
 - 5、轻便、简单、携带方便。
1. Microcomputer intelligent control, freely adjustable operating time.
 2. Eight kinds of automatic mode, three kinds of strength to select.
 3. There are offer two safeguard for your health. The overtemperature protection and the overcurrent protection.
 4. Low voltage design, more security and more comfortable.
 5. Handy, simple, easy to carry.


1. Microcomputer intelligent control, freely adjustable operating time.
2. Eight kinds of automatic mode, three kinds of strength to select.
3. There are offer two safeguard for your health, the over temperature protection and the overcurrent protection.
4. Low voltage design, more security and more comfortable.
5. Handy, simple, easy to carry.

Optional function for all massager bed:

1.3D lifting for upper part of body

4+1 infrared heated jade roller or 6+1 infrared heated jade roller or 8+1 infrared heated jade roller

Jade Massage Bed



2.Thump for leg or back or both

4+1 infrared heated jade roller or 6+1 infrared heated jade roller or 8+1 infrared heated jade roller

Jade Massage Bed



3.Kneading for lower part of body

