

# NLF-200D CPAP System

## Trust Points

- Simplicity: easy to use, easy to move with 5 casters.
- Option: adapt the equipment to your patients and procedures freely.
- Trust: designed and manufactured by Rooe Medical with over 23 years experience in this field.
- Friendly: flexible configurations able to customize your requirements.
- Quality: international standard and advanced technology suitable for all baby patients around the world.
- Convenient: compact interface and big screen give you better operating experience.
- Worldwide: over 4,000 units are installed around the world.

## Features

- Integrated breathing circuit design, ensure easy operating and keep hygiene.
- Built-in air and oxygen blender, ensure stable oxygen concentration.
- Reliable CPAP control valve and pressure monitoring system improve CPAP adjustment precision.

## Safety

- Two casters with brake, ensure stable operating.
- Pressure gauge and flow meter ensure accurate adjusting of parameters.

# NLF-200D CPAP System

*Friendly Powerful Reliable*

## Application


The CPAP System makes a good performance in Neonatal Intensive Care Units (NICU) and other departments. Rooe Medical CPAP System offers you the choices you need in ventilation, monitoring and technique. The expertise in ventilation is based on more than 23 years rich history of developing CPAP solutions that meet your requirements. Great performance is based on high security, high accuracy, high stability and accurate monitoring.



*Cherish your life, Cherish your health!*



Flowmeter



Humidifier



Pressure gauge

## Technical Specifications

Oxygen concentration	Adjusting range: 21% - 100%
CPAP	Adjusting range: 0cmH <sub>2</sub> O - 10cmH <sub>2</sub> O
Continuous flow	Adjusting range: 1.5L/min - 15L/min
Gas source	O <sub>2</sub> , AIR
Pressure	280kPa - 600kPa

Wooden case packing size (main component)	L 700*W 710*H 380mm
G.W.	35KG
N.W.	17.7KG
Wooden case packing size (Trolley)	L 160*W 195*H 1100mm
G.W.	6.7KG
N.W.	2KG

江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19,East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663      Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com      Http://wanrooemed.com