

### **BSS200 Stainless Steel Dressing Trolley**

#### **Benefits**

- Lower Shelf 300mm from the floor to support OH&S requirements to minimize bending/lifting exposure
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier

#### **Specifications:**

Part Number	BSS200
Item	Dressing Trolley-without Drawer
Shelf	2 Shelves, 4*rails on both shelf
Castor	4*100mm, two with brakes
Construction	Stainless steel/Welded
Measurement(mm)	Option 1: 850W*500D*900H Option 2: 700W*450D*900H Option 3: 600W*400D*900H Or other customized size



### **BSS200D/200DB/201D Stainless Steel Dressing Trolley with one Drawer**

		
BSS200D	BSS200DB	BSS201D

#### **Benefits**

- Lower Shelf 300mm from the floor to support OH&S requirements to minimize bending/lifting exposure
- Full Extension Ball Bearing drawer runners which are designed to protect contents impact damage
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Drawers will not fly open when carts are being moved around
- Handles rounded and entirely smooth to minimize infection control and make grasping easier(BSS201D)

### Specifications:

Part Number	BSS200D	BSS200DB	BSS201D
Item	Dressing Trolley-one Drawer Dressing Trolley-one Drawer		
Shelf	2 Shelves, 3*rails on both shelf	2 Shelves, 3*rails on both shelf	2 Shelves
Bucket	none	1 pc	none
Castor	4*100mm, two with brakes 4*100mm, two with brakes		
Construction	Stainless steel/Welded		
Measurement(mm)	Option 1: 700W*450D*900H Option 2: 600W*400D*900H Or other customized size	Option 1: 700W*450D*900H Option 2: 600W*400D*900H Or other customized size	600*360*850

### BSS200D2B Stainless Steel Dressing Trolley with two Drawer

#### Benefits

- Lower Shelf 300mm from the floor to support OH&S requirements to minimize bending/lifting exposure
- Full Extension Ball Bearing drawer runners which are designed to protect contents impact damage
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Drawers will not fly open when carts are being moved around

### Specifications:

Part Number	BSS200D2B
Item	Dressing Trolley-double Drawer
Bucket	1 pc
Shelf	2 Shelves, 4*rails on both shelf
Castor	4*100mm, two with brakes
Construction	Stainless steel/Welded
Measurement(mm)	Option 1: 700W*450D*900H Option 2: 600W*400D*900H Or other customized size



### **BSS200D2 Stainless Steel Dressing Trolley with Double Drawer**

#### **Benefits**

- Lower Shelf 300mm from the floor to support OH&S requirements to minimize bending/lifting exposure
- Full Extension Ball Bearing drawer runners which are designed to protect contents impact damage
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Drawers will not fly open when carts are being moved around

#### **Specifications:**

Part Number	BSS200D2
Item	Dressing Trolley-double Drawer
Bucket	1 pc
Shelf	2 Shelves, 4*rails on both shelf
Castor	4*100mm, two with brakes
Construction	Stainless steel/Welded
Measurement(mm)	760*460*860



### **BSS202D2 Spinal Flat Top Trolley**

#### **Benefits**

- Lower Shelf 300mm from the floor to support OH&S requirements to minimize bending lifting exposure
- Full Extension Ball Bearing drawer runners which are designed to protect contents impact damage
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Drawers will not fly open when carts are being moved around
- Solid construction means a longer product life
- Two (2) years warranty

#### **Specifications:**

Part Number	BSS202D2
Item	Dressing Trolley-double Drawer
Shelf	2 Shelves
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	500w x 500d x 900h



### **BSS201DB Medication Cart**



#### **Benefits**

- Full Extension Ball Bearing drawer runners which are designed to protect contents impact damage
- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Drawers will not fly open when carts are being moved around
- Solid construction means a longer product life
- Two (2) years warranty

#### **Specifications:**

Part Number	BSS201DB
Item	Dressing Trolley-one Drawer
Bucket	1 pc
Instrument shelf	1 pc
Shelf	2 Shelves, 3*rails on both shelf
Castor	2 x 100mm Locking, 2 x 100mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	800w x 500d x 1000h

### BSS203 fan-shaped Instrument Table



- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Solid construction means a longer product life
- Two (2) years warranty

#### Specifications:

Part Number	BSS203
Item	Dressing Trolley
Shelf	2 Shelves, 3*rails on both shelf
Castor	5pcs
Construction	Stainless steel/Welded
Measurement(mm)	1370w x 400d x 1000h

### BSS300/301 Three Layer Stainless Steel Instrument Trolley



Benefits

- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- Solid construction means a longer product life
- Two (2) years warranty

**Specifications:**

Part Number	BSS300 / BSS301
Item	Dressing Trolley
Shelf	3 Shelves, 3*rails on both shelf
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	Option 1: 850w x 500d x 900h Option2: 700w x 450d x 900h

**BSS302 Trolley For Infusion Bottles with three shelves**

- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Four Sides with guardrails
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- Solid construction means a longer product life
- Two (2) years warranty

**Specifications:**

Part Number	BSS302
Item	Trolley For Infusion Bottles with three shelves
Shelf	2 Shelves, 4*rails on both shelf
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	850w x 500d x 900h



BSS204 stainless steel diaper trolley



**Specifications:**

Part Number	BSS204
Item	BSS204 stainless steel diaper trolley
Shelf	2 Shelves
Castor	2 x 100mm Locking, 2 x 100mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	550w x 400d x 800h



### BSS020 Stainless Steel Linen Trolley with dust bag

- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- Solid construction means a longer product life
- Two (2) years warranty

#### Specifications:

Part Number	BSS020
Item	<b>BSS020 Stainless Steel Linen Trolley with dust bag</b>
Shelf	3 Shelves on left, right side with dust canvas bag
drawer	1 pc
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm )	940w x 460d x 940h



### BSS021/022/023 Stainless Steel Linen Trolley with dust bag

		
<b>BSS021</b> Dia:588mm, height:910mm made of stainless steel with 3 inch castor,easy to move waterproof cloth,general color in blue	<b>BSS022</b> 470*425*895mm made of stainless steel with 3 inch castor,easy to move waterproof cloth,general color in blue	<b>BSS023</b> Type I:: 1190x560x990mm Type II: 800x560x860mm typeIII: 560x560x860mm made of stainless steel with3 inch castor,easy to move waterproof cloth,general color in blue

**BSS024 Stainless Steel Linen Trolley with dust bag**



- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- with waterproof canvas dust bag
- Solid construction means a longer product life
- Two (2) years warranty

**Specifications:**

Part Number	BSS024
Item	<b>BSS024 Stainless Steel Linen Trolley with dust bag</b>
Shelf	3 Shelves on left, right side with dust canvas bag
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	940w x 460d x 940h



### BSS025/026 Stainless Steel Linen Trolley with double dust bag



BSS025

BSS026

- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- double waterproof dust bag, one in yellow color, one in blue color
- Solid construction means a longer product life
- Two (2) years warranty

#### Specifications:

Part Number	BSS025/026
Item	<b>Stainless Steel Linen Trolley with dust bag</b>
Dust Bag	One yellow color, one blue color
Castor	4 x 100mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	810W x 530D x 900H

### BSS027/028 S Laundry Trolley with Lid



BSS027 Double Laundry Trolley With Lid

BSS028 Single Laundry Trolley With Lid

- stainless steel or steel power coated optional

- with ABS lid, open via foot pedal

Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort

- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- double waterproof dust bag, one in yellow color, one in blue color
- Solid construction means a longer product life
- Two (2) years warranty

#### Specifications:

Part Number	BSS027	BSS028
Item	Double Laundry Trolley with Lid	Single Laundry Trolley with Lid
Dust Bag	One yellow color, one blue color	Blue or yellow color
Castor	4 x 100mm Swivel	4 x 100mm Swivel
Construction	Stainless steel or steel /Welded	Stainless steel or steel /Welded
Measurement(mm)	800W x 590D x 900H	600W*430D*900H

#### **BSS029 Stainless Steel Linen Trolley with dust bag**

- Precision ball bearing castors which enables the weight of the unit to be absorbed by the bearing and then assist the unit being moved with minimum effort
- Handles rounded and entirely smooth to minimize infection control and make grasping easier
- with waterproof canvas dust bag
- Solid construction means a longer product life
- Two (2) years warranty

#### Specifications:

Part Number	BSS029
Item	<b>Stainless Steel Linen Trolley with dust bag</b>
Shelf	4 Shelves on left, right side with dust canvas bag
Castor	2 x 125mm Locking, 2 x 125mm Swivel
Construction	Stainless steel/Welded
Measurement(mm)	1050w x 500d x 1000h



### **BSS030 Cylinder Trolley**



Made of full stainless steel

Size: 850\*500\*900mm

Tubular frame work.

4 castors

with base plate for support the oxygen cylinder

### **BSS031 Double Bowl Stand**



Size:  $\varnothing$  36\*76-126cm

Made of stainless steel

height can adjust from 76-126CM,

double bowl, diameter: 36cm

5pcs 2 inch castors, 2pcs with brake

### **BSS032 Mayo table with bowl**

Size:  $\phi$  36\*76-126cm

Carton packing : 77\*12\*38cm/2pcs

1. G.W. : 10KG
2. The bowl and bowl holder are made of stainless steel;
3. 5 PP legs base, with castor, 2inch, two with brakes
4. height can adjust from 76-126cm;
5. The bowl can use alone for sterilizing;
6. 5pcs 2 inch castors, 2pcs with brake



### **BSS001/002 Mayo trolley**

	
<p><b>BSS001</b> 680*460*940-1400mm 3 inch castor made of stainless steel</p>	<p><b>BSS002</b> 650*400*900-1400mm 3 inch castor made of stainless steel</p>

### **BSS003 Treatment trolley**

made of stainless steel  
overall size: 380\*450\*920mm  
double tray  
double container in plastic  
4 castor, diameter of 100mm



**BSS004/005 Stainless steel kick bucket**



**BSS004**

Inner diameter:260mm  
outside diameter:300mm  
height 240mm  
we can make according to customer's size

**BSS005**

Bucket diameter 35cm, height 43cm

**BSS007 Stainless Steel Bed Pan**




Size: 350\*280\*97  
made of stainless steel  
lid: optional



**BSS006 Stainless Steel waste bin with foot pedal**



根据需要 多种尺寸可供选择



made of stainless steel with foot pedal to open the lid

3/5/8/12/20/30L available

with blue protection film, when in use, please remove it

with lifting handle,easy for move

code	size(CM)	Inner size(CM)	carton(CM)	Qty/carton(PCS)
BSS006-03L	16.7X24.8	21X17X25.5	64X53.5X53.5	18
BSS006-05L	20.5X28.6	22X22X29.8	69X46X61	12
BSS006-08L	22.5X33	25X25X34	78X53X70	12
BSS006-12L	25X39	26.5X26.5X40.5	55.5X55.5X82	8
BSS006-20L	29X44	31X31X44.5	126X64X46	8
BSS006-30L	29X65.5	31X31X66.5	64X64X68.5	4



蓝色的保护膜收到货以后可以剥掉，光亮如镜

特大脚踏设计，使盖开启便捷



桶盖密封设计，能有效的隔绝异味

安装手提吊环，可方便的随意移动