




WANROOE
MEDICAL

www.wanrooemed.com

Height controlled by foot pedal
USA brand Power Packer Oil Pump
Three function/Five Function available



BAH300

Three Function Hydraulic Bed

张家港市万容机械有限公司
WANROOE MACHINERY CO., LTD



WANROOE
MEDICAL

www.wanrooemed.com

BAH300 Three Function Hydraulic Bed



Feature:

1. Bed base, frame and leg are all made of cold-rolled steel plate and tube and coated by electrostatic spray.
2. The quality attains national standard. Head and foot boards are made of imported PP materials.
3. 6-bars siderail
4. with infusion stand, 4 hooks
5. ordinary castor, diameter 5 inch, two with brakes
6. Backrest and leg section gas spring control optional
Trendelenburg and anti-trendelenburg function we can add according to customer requirement.

Material	Steel Framework
Overall size	215*98*45/75cm Patient Touching Surface 90*196CM
Hi-Low	45-75cm by foot pedal
Backrest lifting angle	0-75°
Knee part lifting angle	0-40°
Silent and Dustproof Castor diameter	125mm
loading	250kgs



Tel:0086-512-58810860 Fax:0086-512-58810663
Email: info@wanrooemed.com www.wanrooemed.com