

BDT001A/BDT001B/BDT001C/BDT001D Over Bed Table, steel column



A

B

C



D

E

F

Features:

1. Tabletop: A type made of plastic PP blow molding,
B type: Hospital movable over bed table with MDF surface
C type: plastic PP blow molding with bottle holder on tabletop
D type: plastic PP blow molding with bottle holder on tabletop, with drawer
2. Lifting Column: steel with plastic coating
3. Base made of steel with plastic coating
4. Four castors, two with brakes, easy to move, durable

Specifications:

Item code	Overall size	G.W/N.W(kg)	Qty/package	Package size
BDT001A	900*450*750-1020mm	15/13.5	1	900*560*170mm
BDT001B	900*450*750-1020mm	17.5/16.5	1	900*560*170mm
BDT001C	900*450*750-1020mm	16/14.5	1	900*560*170mm
BDT001D	900*450*750-1020mm	16.5/15	1	900*560*170mm
BDT001E	775*385*775-1055mm	16.5/15	1	900*560*170mm
BDT001F	900*450*750-1020mm	20/18	1	900*560*200mm

BDT001G OVER BED TABLE



1. ABS material table board, no toxic effect on body & be easy to clean, general color: **white or blue**
2. Two sections ABS tabletop, adjustable angle for reading books etc.
3. Steel column
4. The height can adjustable: 74cm-102cm
5. With 4pcs castor for move easy, two with brakes
6. Press knock switch, adjust height by your demanded.

Dimensions:

Length: 900*450*750-1020mm

BDT002A/BDT002B Over Bed Table, aluminum alloy lifting column



BDT002A

BDT002B

Features:

1. Tabletop: A type made of plastic PP blow molding, B type: Hospital movable over bed table with MDF surface, aluminum packed edge.
2. Lifting column: aluminum profile with electrostatic spraying
3. Base made of steel with plastic coating
4. Four castors, two with brakes, easy to move, durable

WANROOEMED reserves the right to change specifications without notice.

Specifications:

Item code	Overall size	G.W/N.W(kg)	Qty/package	Package size
BDT002A	890*440*830~1040MM	22.7/21.7	2	915*465*320mm
BDT002B	780*386*830~1040MM	22.7/21.7	2	820*420*200 mm

BDT002C Over Bed Table



1. ABS material table board, non-toxic, easy to clean.
2. Size: 900*450*930-1230mm
3. Steel column
4. The height can adjustable: 93cm-123cm
5. 4 castors for move easy
6. Press knock switch, adjust height by your demand

BDT003A/BDT003B Stainless Steel Over Bed Table



Specifications:

Item code	Overall size	G.W/N.W(kg)	Qty/package	Package size
BDT003A	800*400*780-1020mm	22.7/21.7	1	920x475x120mm
BDT003B	900*450*700-1000mm	13/12	1	800x460x120mm